

Phase 2 – Maintenance Phase

As soon as the extremities have been decongested and the circumference fluctuations are brought under control from MLD therapy, medical care experts advise patients to wear medical compression garments to manage the condition. Juzo has a full line of medical compression garments to help manage every stage of lymphedema.



Juzo Soft & Dynamic Compression Garments

Comfortable Circular Knit without Seam

The compressive effect of Juzo circular knit helps manage the reduced circumference achieved by MLD. Circular knit garments help preserve the size of the limb for optimal lymphedema management. Juzo's circular knit products are extremely comfortable, durable & great for mild to moderate swelling.

- Medical containment
- Not made from natural rubber
- Machine wash and dry
- Comfortable and breathable
- Seamless
- Many options, styles & sizes



Juzo Expert & Strong Compression Garments

Comfortable Flat Knit with Seam

Flat knitted custom compression garments help ensure the long-term success of decongestive therapy. Flat knit garments have a high degree of containment and are customized to any shape or size. Juzo's flat knit products are extremely comfortable and durable. The flat knitted garments offer maximum therapeutic results for moderate to extreme swelling.

- Maximum containment
- Not made from natural rubber
- Machine wash and dry
- Flat and non-irritating seam
- Any size and any shape

Juzo Compression Solutions

For Lymphedema



Look for Juzo. Find Comfort.

Juzo USA
P.O. Box 1088
Cuyahoga Falls, OH 44223
www.juzousa.com

Juzo Canada Ltd.
1100 Burloak Dr. Ste. 300
Burlington, ON L7L 6B2
www.juzo.ca

© Juzo is a registered trademark of Julius Zorn, Inc. © 2017 Julius Zorn, Inc.

BRLYMPHTREAT-0917



Freedom in Motion




Freedom in Motion

Effective Management of Lymphedema

Lymphedema is a condition that occurs when there is an accumulation of lymphatic fluid in the tissue caused when lymph vessels are unable to adequately drain lymph fluid. The international gold standard treatment option for lymphedema is Complete Decongestive Therapy (CDT), which involves manual lymphatic drainage (MLD), compression management, consistent skin care and exercise. CDT occurs in two phases: phase one, a decongestion phase and phase two, a maintenance phase.

Juzo offers high-quality medical compression garments that support patients during both phases of CDT.



Phase 1 – Decongestion Phase

During phase one, a lymphedema therapist performs MLD on the patient so accumulated lymph fluid can drain through proper channels. The patient is then advised to wear lymphedema compression bandages. An effective supplement to this treatment is the use of short-stretch compression wraps. They can be adjusted for the individual patient and are suitable as a reduction tool in the management of swelling changes during the phase one MLD process.



Juzo SoftCompress Bandage & Gauze

Compression Aid

The short stretch bandage consists of 100% cotton that is comfortable on the skin and machine washable. Juzo short stretch bandage exerts low resting and high working pressure. The gauze is great for sensitive skin, joints and hard to bandage areas such as genitalia, fingers and toes.

- Machine wash
- Bandage 100% cotton and skin friendly
- Gauze high elasticity and easy to use



Juzo SoftCompress

Compression Bandage Liner

The comfortable Juzo SoftCompress bandaging liner and pads can be worn under short-stretch bandages or wraps and have a skin-friendly cotton surface. The vertically stitched foam channels create differences in pressure, thereby promoting lymph drainage. Juzo SoftCompress allows for a quicker and more efficient application of bandages or wraps.

- 100% cotton surface
- Fast and easy
- Improved outcomes



Juzo Therapy Tape™

Maximum Performance

Kinesiology tape helps the skin to support and stabilize certain muscles and joints without restricting the body's range of motion. When the tape is applied it assists with stimulating the movement of lymph by stretching the skin wherever it is placed and more effectively increasing the flow of lymphatic fluid.

- Adhesive that can be worn in any condition or environment including water, for up to 7 days
- Can be worn 24/7 to help you stay active & healthy
- Can be used as an adjunct therapy



Juzo Compression Wrap

Short Stretch Compression Wrap

The Juzo Compression Wrap is a short stretch garment that can be employed as an alternative to conventional multi-layer bandages. Thanks to its short stretch properties, it delivers high working pressure which will help evacuate fluid out of the limb. The Juzo Compression Wrap may also be used as a compliment to medical compression garments. The fasteners make it easy to apply, and easy to adjust for fluctuations in body circumference. The special notches and minimal overlying straps resist gapping and bulk. The wrap is available in foot, calf, knee, thigh, arm and hand styles.

- Reversible color, like getting two wraps in one
- Short stretch for high working pressure
- Pull tab aid for hassle free application

New Styles: Knee, Thigh, Arm & Hand

Juzo Accessories for Decongestive Phase

Pads and foam targeted for an increase in pressure

Juzo SoftCompress pads and Juzo TheraPad™ are effective aids for lymphedema. They are especially suitable for the targeted application of increased pressure and gently stimulating the tissue, for specific area of edema to reduce fibrosis and encourage lymph drainage.

Juzo SoftCompress Pads



JuzoTheraPad™

