

PRODUCT RANGE

Phlebology • Lymphology • Scar therapy



Our Customer Service department – always there for you

Customer service is a top priority at Juzo – and has been for more than 100 years. Our passionate, reliable and flexible staff is happy to advise you on any aspect of compression garments. Your contact persons are knowledgeable in all areas of our quality products and are glad to assist you with any questions or issues you might have.

Our Customer Service department offers ...

- general information
- individual product advice
- order acceptance
- feedback on the status of your order
- acceptance of alteration or repair orders

Our Technical Service department can advise you on ...

- specialist questions regarding flat knitting and oedema treatment
- complex and complicated flat knitted garments



Edith Männer
Head of Customer Service



Markus Daferner
Head of Technical Service

You can contact us by phone
Monday to Friday.

Please feel free to send us
an e-mail or a fax.

Juzo Customer Service International

Office hours 7:00 am – 6:00 pm
Phone +49 (0)8251 901 300
Fax +49 (0)8251 901 321
E-mail order@juzo.com

Juzo Customer Service UK

Office hours 9:00 am – 5:00 pm
Phone +44 (0)1382 82 66 20
Fax +44 (0)1382 82 66 41
E-mail sales@juzo.co.uk

More than 100 years of experience in compression therapy

The history of the family-run business Juzo extends back to the year 1912. It was in that year that Julius Zorn GmbH (Juzo) was founded in Zeulenroda. This small city in the southeast of Thuringia was already regarded at the time as the elastic stocking manufacturing capital – a place where a high level of know-how and experienced employees came together. The hosiery industry in Zeulenroda has its roots as far back as 1750 and gained an international

reputation in the 19th century. In the year 1919 Juzo revolutionized its manufacturing processes: for the first time ever the stockings were now knitted rather than woven. The demand for these novel products increased steadily thereafter. By the 30s Juzo was already exporting its compression products to many countries in the world. After the Second World War, the grandson of the founder Julius Zorn settled in Aichach and developed the company

Juzo into an international business. Juzo introduced seamless compression stockings in 1951, and in 1962 Juzo began manufacturing the first latex-free compression stockings. Thanks to its consistent focus on quality and further developing its products, Juzo has been successfully growing its international business and currently employs around 900 people all over the world. The manufacturing facilities at the company's site in Aichach are the most modern of their kind.

Summary



Compression therapy

| | |
|-------------------------|----|
| Versatile fabrics | 04 |
| Product characteristics | 06 |
| Quality standards | 11 |




Phlebology

| | | | |
|---------------------------------|----|------------------------------|----|
| The fine fashionable stockings | 12 | The experts with seam | 20 |
| The sporty allround stockings | 14 | Varodem® Hospital | 42 |
| The intensive therapy stockings | 16 | Juzo® Compression Wrap | 43 |
| Juzo® Helasana® | 18 | Juzo® Liner Stockings Silver | 52 |
| Juzo® Ulcer Pro | 19 | | |



Lymphology

| | | | |
|---------------------------------|----|--------------------------------------|----|
| The intensive therapy stockings | 16 | Compression gloves | 34 |
| The experts with seam | 20 | Compression segments thorax | 38 |
| The strong experts with seam | 22 | Compression segments head | 40 |
| Compression segments leg | 24 | Juzo® ACS light | 44 |
| Compression foot-toe-portions | 24 | Juzo® SoftCompress | 46 |
| Compression stump shrinkers | 26 | Juzo® Postoperative Compression Bras | 50 |
| Compression sleeves | 30 | Juzo® Liner Stockings Silver | 52 |



Scar therapy

| | | | |
|-------------------------------|----|--------------------------------------|----|
| The experts with seam | 20 | Compression segments thorax | 38 |
| Compression segments leg | 24 | Compression segments head | 40 |
| Compression foot-toe-portions | 24 | Juzo® Postoperative Compression Bras | 50 |
| Compression stump shrinkers | 26 | Juzo® ScarPad silicone patch | 51 |
| Compression sleeves | 30 | Juzo® Liner Stockings Silver | 52 |
| Compression gloves | 34 | | |



Special products

| | | | |
|--------------------------------------|----|------------------------------|----|
| Varodem® Hospital | 42 | Juzo® ScarPad silicone patch | 51 |
| Juzo® Compression Wrap | 43 | Juzo® Liner Stockings Silver | 52 |
| Juzo® ACS light | 44 | Juzo® Special Socks | 53 |
| Juzo® SoftCompress | 46 | Juzo® Light Line | 54 |
| Juzo® Postoperative Compression Bras | 50 | | |



Accessories

| | | | |
|---|----|---------------------|----|
| Juzo® Slip-on and Taking-off Aids | 56 | Juzo® Pressure Pads | 60 |
| Juzo® Arion Slip-on and Taking-off Aids | 57 | Juzo® Lymphpads | 60 |
| Juzo® Care and Cleaning Products | 59 | | |



Versatile fabrics For effective compression therapy

Medical compression stockings from Juzo guarantee correct pressure values and a well-defined pressure gradient. The pressure decreases continuously from the ankle section upwards towards the heart, and also towards the toes. We can provide a wide variety of garments, from below-knee stockings to pantyhose, from glove to face mask, depending on the therapeutic application.

Compression fabrics are categorized into 4 compression classes (CCL) based on the level of pressure.

Possible indications can be, for example:

18 – 21 mmHg (CCL 1) // mild compression

- minor obstructive disorders accompanied by feelings of tiredness or heaviness in the legs
- isolated varicose vein formation
- preventative measure during pregnancy
- therapeutic treatment of scars

23 – 32 mmHg (CCL 2) // moderate compression

- varicose veins
- leg swelling during pregnancy
- after varicose vein sclerosis or operations
- for the long-term treatment of healed leg ulcers
- propensity for oedema formation exists

34 – 46 mmHg (CCL 3) // strong compression

- varicose veins accompanied by damage to the deep leg veins
- after healing of venous crural ulcers
- when a severe swelling tendency exists
- chronic venous insufficiency
- severe propensity for oedema formation

≥ 49 mmHg (CCL 4) // very strong compression

- severe lymphatic obstruction
- elephantiasis
- lymphoedemas with fibrosis



Use in phlebology

Circular knitted seamless products are predominantly used in venous therapies. Today's seamless compression stockings are often indistinguishable from normal fine-mesh stockings and comfortable to wear. Of course, they still have the necessary compression properties to guarantee an effective therapy. These compression stockings are used for indications such as varicosity, chronic venous insufficiency, or post-operatively after varicose vein stripping.

How it works

The compressive fabric exerts a pressure on the leg, and thereby also on the veins. This reduces the diameter of the veins, thus enabling the damaged venous valves to close again. The blood flows out of the legs better, and painful swelling and obstructions recede.

Circular knitted products from Juzo: Juzo Attractive, Juzo Hostess, Juzo Soft, Juzo Soft Rib, Juzo Dynamic Cotton Rib, Juzo Dynamic, Juzo Dynamic Silver, Juzo Dynamic Cotton and Juzo Helasana



New knitting technology

An innovative seamless knitting technology enables the seamless manufacture of compression garments. It eliminates seams completely. The soft 3-dimensional fabric conforms perfectly to anatomical features and produces a well-defined pressure gradient.

Introduction to oedema therapy

Seamless knitted products, for example the Juzo Classic Seamless, are especially suitable for treating incipient lymphoedema where the anatomical situation is still normal.

Seamless knitted products from Juzo: Juzo Classic Seamless



Use in lymphology

Flat knitted compression products are the ideal complement to lymph drainage treatment when treating oedemas. Compression garments help secure the long-term success of decongestive therapy and they are generally indispensable in the area of lymphology. The thicker seamed fabrics, e.g. Juzo Expert Strong, offer additional benefits, for example a massage effect that promotes lymph drainage. Furthermore, these fabrics do not yield to the oedema. In conjunction with movement, they produce a high therapeutic pressure that provides optimum compression of the oedematous tissue.

Use in scar therapy

Compression garments can also contribute to very good treatment results in the area of scar therapy. To ensure no proliferative scarring occurs, it is essential to apply external compression (18 – 32 mmHg, CCL 1 – 2) as soon as the wound has closed. This accelerates scar maturation and reduces the thickness of the scar.

Custom-made fabrics

The flat knitted products Juzo Expert and Juzo Expert Strong are knitted stitch by stitch to precisely match the body measurements of the patient, and come in various styles to suit the particular therapeutic application. Juzo garments are made from different combinations of materials, for example incorporating cotton or silver, and are therefore able to meet a wide variety of patient requirements.

Pure silver with antibacterial effect

Patients affected by burns, scalding or scarring are especially prone to infections by pyogenic bacteria (e.g. *Staphylococcus aureus*). Our special-purpose compression fabric with a continuous knitted-in silver thread (e.g. Juzo Expert Silver) has a natural antibacterial effect on these bacteria on the skin. Silver-containing fabrics can therefore restore the natural balance of inflamed or irritated skin as well as reduce itchiness, thereby promoting better healing.

Custom-made garments with cotton

Compression products with an especially high cotton content are extremely soft and very smooth because the skin is in contact with pure cotton. Thanks to this softness, fabrics like Juzo Expert Cotton are especially recommended for delicate and sensitive skin.

Flat knitted products from Juzo: Juzo Expert, Juzo Expert Silver, Juzo Expert Cotton, Juzo Expert Strong, Juzo Expert Strong Silver and Juzo Expert Strong Cotton

Product characteristics

Circular knitting (seamless)



Perfect elasticity

The exceptional stretch properties guarantee an outstanding fit and high wearing comfort.



For all styles

A large colour selection with additional trend colours, knitted-in patterns, or batik look – there is something to suit every style.



Versatile silicone border options

An optimum hold for every requirement – whether for the fashion-conscious woman, an active man or for sensitive skin.



Top border

Below-knee stockings and sleeves with smooth, non-constricting and comfortable-to-wear top border.



Variety of standard sizes

Almost all models are available in standard sizes – for special needs custom-made is possible.



Flat knitting (with seam)



Flat knitting



Versatile fabrics

Special styles offering a high level of skin compatibility, also with silver thread or cotton.



Colour variety

Seven fabulous standard colours that are outstandingly combinable – changing trend colours and a batik look round off the colour range.



Flat seams

Comfortable to wear thanks to the flat, elastic seam.



Extra options

Angled openings, pressure pads, instep stitching and lots more extra options guarantee a perfect fit.



Full-body garments

Custom-made items from toe parts to face masks with optimum fit.



Flat knitting (with seam on the outside)



Soft fabrics

Optimum sun protection thanks to its certified UV protection factor of 80 – special styles also with silver.




Colour variety

Seven fabulous standard colours that are outstandingly combinable – changing trend colours and a batik look round off the colour range.



Seams on the outside

The flat, elastic seams on the outside are barely noticeable – a very soft transition for the delicate healing skin.



Extra options

An inelastic back zone, angled openings, prosthesis pockets, instep stitching and other extra options for even more comfort.



Full-body garments

Custom-made items from toe parts to face masks with optimum fit.



Iron-on patches for children




Especially for custom-made children's garments

Our iron-on patches especially for all custom-made children's leg, arm, hand, thorax and head garments in Juzo Expert and Juzo Expert Strong (all colours) as well as Juzo Expert Silver, Juzo Expert Cotton, Juzo Expert Strong Silver and Juzo Expert Strong Cotton (colour Almond).

Here you can find a selection of iron-on patches for the children's compression products of "Julius" and his friends.


You can find all patches on www.juzo.com/iron-on-patches



"Lollipop"



"Diving suit"



"Whale"



"Chair"



"Henriette"



"Fridolin"



"Ferdinand"



"Gustav"



"Lotte"




"Gottfried"



"Seagrassflipper"

Silicone borders



Top border

Top border without silicone – like the usual top of a sock.
For fastening of thigh stockings (without silicon border) adhesive lotion or another attachment possibility may be necessary.



Elastic border

Durable border without silicone beads (dots).
It stabilises the top border but has no adhesive properties and used for AD stockings and sleeves.



Silicone border

The silicone beads (dots) provide additional skin-friendly adhesion and prevent the stocking from slipping. The small beads allow the air to circulate well. The discreet silicone border is also suitable for male compression stocking wearers.

For AD, CG and SG: thin silicone border (3.5 cm)



Decorative silicone border

The decorative silicone border has the same characteristics as the silicone border. This border option is especially demanded by women because there is a pretty but inconspicuous pattern on the outside of the border.



Lace border

A border option with a 3-row, continuous circular silicone coating on the inside. The silicone has excellent skin-friendly adhesion characteristics. Ladies appreciate the fine-socking look and the elegant appearance due to the woven in floral pattern.



Comfort silicone border

A unique, sensitive silicone coating provides a comfortable fit and is very gentle on the skin. The border is very flexible and fits perfectly on the thigh. The breathable surface reduces perspiration. The comfort silicone border has a low profile – especially designed for fashion-conscious and ambitious women.

Quality standards

Verified quality



Oeko-Tex® Standard 100

“Oeko-Tex Standard 100” is a globally-standardized and independent testing and certification system that applies to all stages of processing, from raw materials through to the intermediate and final product. The more intensive the skin contact, the stricter the requirements that must be met.

All Juzo compression fabrics, both the circular knitted and flat knitted varieties, have been awarded the “Oeko-Tex Standard 100” label for textiles. This enables the company Juzo to guarantee that all components that go into a product are completely harmless to humans and the environment. The overall requirements of this standard far exceed those of existing national legislation.



Quality seal RAL

All medical compression stockings from Juzo comply with the relevant RAL standard. The test specifications contain a comprehensive set of binding guidelines on the manufacture of medical compression stockings and sleeves, for example:

- material characteristics (harmless to humans)
- periodic quality inspections, for example by the Hohenstein Research Institute
- rubber-elastic properties in lengthwise and crosswise directions
- compression classes 1 – 4, for sleeves compression classes 1 – 3
- manufacture of custom-made and standard sized stockings



UV Standard 801

Juzo Expert in 18 – 32 mmHg (CCL 1 – 2) offers optimum sun protection and a UV protection factor of 80 – as verified and certified by the Hohenstein Institute.



ISO certification

Juzo implemented its general quality management system back in 1998, and was one of the first companies in this industry to be certified to EN ISO 9001 and EN 46001. Juzo is currently certified to EN ISO 9001 and EN ISO 13485, which superseded EN 46001.

The fine fashionable stockings



Juzo® Attractive


Glossy with a fine-stocking look

The velvety transparency and exceptionally fine knit make this compression stocking a real trendsetter. Thanks to its fine-stocking look, Juzo Attractive is perfectly suited for everyday use and is a versatile addition to your wardrobe. The fashion highlight of this range is the "Vision" pattern, which can be incorporated on request. The compression pantyhose with a special pregnancy body part optimally supports both the expectant mother and the baby.

"Vision" pattern
available in standard sizes
with closed toe

Closed toes
in standard sizes suitable for
shoe sizes 37 – 42 (UK: 4 – 8)

Open toes
in standard sizes suitable for
shoe sizes 37 – 42 (UK: 4 – 8)



Juzo® Hostess®

Airy feeling with transparent look

Thanks to its unique weave, Juzo Hostess is one of the most breathable products in the fine compression stockings segment. The fibres are treated in a special process that makes the stocking especially permeable to air as well as highly durable.

Closed toes
in standard sizes suitable for
shoe sizes 37 – 42 (UK: 4 – 8)

Open toes
in standard sizes suitable for
shoe sizes 37 – 42 (UK: 4 – 8)



Juzo® Soft

Super soft for feel good moments

The fine microfibre fabric with its optimum moisture transport guarantees a comfortable skin environment and a pleasant, soft feeling on the skin. Juzo Soft combines refined fashion with the ultimate in wearing comfort. You can make a fashion statement with the Juzo Soft Leggings which can be combined with every outfit – just feel your best!

Closed toes
in standard sizes suitable for
shoe sizes 37 – 42 (UK: 4 – 8)

Open toes
in standard sizes suitable for
shoe sizes 37 – 42 (UK: 4 – 8)




| Model | Juzo® Attractive | | Juzo® Hostess® | | Juzo® Soft |
|---|---|--------------|----------------|--------------|-----------------------------|
| | Article number / Compression class | 2701 2702 | 2781 2782 | 2501 2502 | 2581 2582 |
| AD normal / short / extra short | | ● | | ● | ● |
| | AD normal / short / extra short with thin silicone border | ● | | ● | |
| AF | AF | | | | ● |
| | AF with silicone border | ● | | ● | ● |
| | AF with decorative silicone border | ● | | ● | ● |
| | AF with lace border | ● | | ● | ● |
| | AF with comfort silicone border | ● | | ● | ● |
| AG normal / short / extra short | AG normal / short / extra short | | | | ● |
| | AG normal / short / extra short with silicone border | ● | | ● | ● |
| | AG normal / short / extra short with wide silicone border | ● | | ● | ● |
| | AG normal / short / extra short with decorative silicone border | ● | | ● | ● |
| | AG normal / short / extra short with wide decorative silicone border | ● | | ● | ● |
| | AG normal / short / extra short with lace border | ● | | ● | ● |
| | AG normal / short / extra short with comfort silicone border | ● | | ● | ● |
| | AG normal / short with hip attachment | | | | ● |
| AT normal / short / extra short with elastic body part | AT normal / short / extra short with elastic body part | ● | | | |
| | AT normal / short / extra short with highly-elastic body part | | | | ● |
| | ATU normal / short / extra short with pregnancy body part | | ● | | |
| | AT with open crotch | | | | ● |
| | AT with fly | | | | ● |
| BT normal / short with highly-elastic body part | | | | | ● CCL 1, not in custom-made |

Standard colours


● available in standard sizes I – VI and custom-made ● only available in custom-made

Juzo Attractive, Juzo Hostess and Juzo Soft



Trend colours · Batik Collection

Juzo Attractive and Juzo Soft



See the "Product overview" for the current trend colours and Batik Collection. The Batik Collection is available in all trend and standard colours (except Almond) in tie dye in combination with white or black.

The sporty allround stockings



soft fabric



ribbed structure

Juzo® Soft Rib

Comfortably soft with a ribbed structure

The striking ribbed structure is what gives Juzo Soft Rib its unique look. Thanks to its outstanding wearing characteristics, this stocking is ideal for everyday wearing – both for men and women alike. A tangibly soft and gentle wearing sensation.



Closed toes

in standard sizes suitable for shoe sizes 38 – 44 (UK: 5 – 9.5)



sport



ribbed structure



with cotton

Juzo® Dynamic Cotton Rib

The allrounder for active people

This stocking is perfectly tailored for the requirements of active people. With its attractive look, stable and durable structure, and opaque fabric, this stocking has been designed for use in a wide range of outdoor activities and sports. Thanks to the high cotton content of Juzo Dynamic Cotton Rib, it is also exceptionally skin friendly, making it comfortable to wear in the office as well.



Closed toes

in standard sizes suitable for shoe sizes 39 – 45 (UK: 5.5 – 10)
short foot part for shoe shoe sizes 36 – 39 (UK: 3.5 – 5.5)



| Model | Juzo® Soft Rib | Juzo® Dynamic Cotton Rib |
|---|----------------|--------------------------|
| Article number / Compression class | 2021 2022 | 3521 3522 |
|  AD normal / short / extra short | ● | ● |
| AD short / extra short with short foot part | | ● |
| AD normal / short / extra short with thin silicone border | ● | ● |
| AD short / extra short with short foot part and thin silicone border | | ● |

● available in standard sizes I – VI and custom-made

Standard colours

Juzo Soft Rib



Juzo Dynamic Cotton Rib



The intensive therapy stockings



 durable

Juzo® Dynamic

Durability you can count on

Juzo Dynamic, with its tight and unobtrusive weave, is visually discreet. It offers a perfect fit, and retains its form even during high levels of physical exertion. This makes it the ideal circular knitted stocking for demanding requirements.



Closed toes
in standard sizes suitable for shoe sizes 37 – 42 (UK: 4 – 8)



Open toes
in standard sizes suitable for shoe sizes 37 – 42 (UK: 4 – 8)



 with silver

 durable

Juzo® Dynamic Silver

With antibacterial action to combat skin irritations

The continuous knitted-in silver thread with antibacterial effect assists the compression therapy by natural means. This compression stocking restores the natural balance of inflamed or irritated skin and noticeably reduces itchiness.



Closed toes
only in custom-made



Open toes
only in custom-made



 with cotton

 durable

Juzo® Dynamic Cotton

With cotton for skin friendliness

This stocking, with comfortable cotton on the inside, is especially skin-friendly and feels pleasant on the skin even with skin irritations and allergies. The cotton, a natural fibre, optimally supports the functional characteristics of Juzo Dynamic.



Closed toes
in standard sizes suitable for shoe sizes 37 – 42 (UK: 4 – 8)



Open toes
in standard sizes suitable for shoe sizes 37 – 42 (UK: 4 – 8)




| Model | Juzo® Dynamic | Juzo® Dynamic Silver | Juzo® Dynamic Cotton |
|---|----------------------|----------------------|----------------------|
| Article number / Compression class | 3511 3512 3513 | 3511 3512 3513 | 3511 3512 |
|  AD normal / short / extra short | ● | ● | ● |
| AD normal / short / extra short with thin silicone border | ● | ● | ● |
|  AF | ● | ● | ● |
| AF with silicone border | ● | ● | ● |
| AF with decorative silicone border | ● | ● | ● |
| AF with lace border | ● | ● | ● |
|  AG normal / short / extra short | ● | ● | ● |
| AG normal / short / extra short with silicone border | ● | ● | ● |
| AG normal / short / extra short with wide silicone border | ● | ● | ● |
| AG normal / short / extra short with decorative silicone border | ● | ● | ● |
| AG normal / short / extra short with wide decorative silicone border | ● | ● | ● |
| AG normal / short / extra short with lace border | ● | ● | ● |
| AG normal / short with hip attachment | ● | ● | ● |
| AG normal / short with body part (Expert) | ● | ● | ● |
|  AT normal / short / extra short with light compressive body part | ● | ● | ● |
| AT one-legged pantyhose with leg extension | ● | ● | ● |
| AT with open crotch | ● | ● | ● |
| AT with fly | ● | ● | ● |
| AT with body part like knit of stocking | ● CCL 2 – 3 | ● CCL 2 – 3 | ● CCL 2 |

● available in standard sizes I – VI and custom-made ● only available in custom-made

Standard colours

Juzo Dynamic



Juzo Dynamic Silver



Juzo Dynamic Cotton





Juzo® Helasana®

The cost-effective alternative

The supple and elastic stocking offers value for money. Despite its fine knit it is durable. The stocking is very breathable. It is made of helanca and a one-way stretch compression stocking.



Open toes


in standard sizes suitable for shoe sizes 37 – 41 (UK: 4 – 7.5)

| Model | Juzo® Helasana® |
|---|-----------------|
| Article number / Compression class | 4512 |
|  AD normal / short | ● |
|  AG normal / short | ● |
|  AG normal / short with silicone border | ● |
|  AG normal / short with hip attachment | ● |

● available in standard sizes II – V

Standard colour





Juzo® Ulcer Pro
For the treatment of crural ulcers

Using its long-standing experience as a manufacturer of compression garments and working in close collaboration with medical technologists and therapists, Juzo has developed a compression stocking system that makes it easier to treat patients with ulcus cruris venosum.

The treatment set comprises two types of compression stockings that are used in combination – one Ulcer Expert and two Ulcer Liners which hold wound dressings and pads securely in place.

The compression stocking system has a lot of compelling features, amongst others it fits securely even during exercise and it can be worn with normal shoes.



Juzo Ulcer Liner



Juzo Ulcer Expert

Juzo® Ulcer Liner

The circular knitted Juzo Ulcer Liner stocking is suitable for long-term compression – 24 hours a day.




Closed toes
in standard sizes suitable for shoe sizes 37 – 42 (UK: 4 – 8)


Juzo® Ulcer Expert

The flat knitted Juzo Ulcer Expert cover stocking, enhances the compression effect during the active phase. It is worn over the liner stocking during the day to increase the level of compression.




Open toes
in standard sizes suitable for shoe sizes 37 – 42 (UK: 4 – 8)

| Model | Juzo® Ulcer Pro |
|---|-----------------|
| Article number / Compression class | 3073 |
|  Compression stocking system AD normal / short set with 2 Ulcer Liner and 1 Ulcer Expert | ● |

| Model | Juzo® Ulcer Pro Add |
|---|---------------------|
| Article number / Compression class | 2571 |
|  Replacement set AD normal / short standard size: set with 3 Ulcer Liner custom-made: set with 2 Ulcer Liner | ● |

Standard colours



● available in standard sizes 1 – 6 and in custom-made

For further information regarding the Juzo Ulcer Pro please see the brochure "Ulcus cruris venosum".

The experts with seam



 UV protection

 flat seam



Juzo® Expert

Pleasantly soft for tailored treatment

Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement and is very comfortable to wear. This special fabric is ideally suited for full-body garments and, thanks to its certified UV protection factor of 80 (18 – 32 mmHg, CCL 1 – 2), is also ideal for use in scar therapy.



Closed toes
only in custom-made



Open toes
in standard sizes suitable for
shoe sizes 37 – 41 (UK: 4 – 7.5)
long foot part for shoe sizes 40 – 44 (UK: 7 – 9.5)



 with silver

 flat seam



Juzo® Expert Silver

Soft fabric with antibacterial effect


Thanks to the knitted-in silver thread, this fabric provides not only compression but also a natural antibacterial effect. This restores the natural balance of inflamed or irritated skin, thereby reducing itchiness.



Closed toes
only in custom-made



Open toes
only in custom-made



 with cotton

 flat seam



Juzo® Expert Cotton

Natural fibres for extra softness

This stocking is made from almost 50 percent of cotton and covers the skin with ultrapure cotton. The natural fibres make it extremely soft and smooth. Juzo Expert Cotton is therefore especially comfortable to wear for sensitive, delicate skin.



Closed toes
only in custom-made



Open toes
in standard sizes suitable for
shoe sizes 37 – 41 (UK: 4 – 7.5)
long foot part for shoe sizes 40 – 44 (UK: 7 – 9.5)



| Model | Juzo® Expert | Juzo® Expert Silver | Juzo® Expert Cotton |
|--|-------------------------------------|------------------------------|----------------------|
| Article number / Compression class | 3021 3022 3023 3024 | 3021 3022 3023 3024 | 3021 3022 3023 |
| AD normal / short | ● CCL 1 – 2 AD; ● CCL 2 AD short | ● | ● CCL 1 – 2 AD |
| AD normal max / short max | ● CCL 2, not in custom-made | | |
| AD normal / short with thin silicone border | ● CCL 2 | ● | ● CCL 2 AD |
| AD normal max / short max with thin silicone border | ● CCL 2, not in custom-made | | |
| AD normal / short with elastic border | ● CCL 2 | ● | ● CCL 2 AD |
| AF | ● | ● | ● |
| AF with silicone border | ● | ● | ● |
| AF with decorative silicone border | ● | ● | ● |
| AF with lace border | ● | ● | ● |
| AG | ● | ● | ● |
| AG with silicone border | ● | ● | ● |
| AG with decorative silicone border | ● | ● | ● |
| AG with lace border | ● | ● | ● |
| AG with hip attachment | ● | ● | ● |
| AT | ● | ● | ● |
| AT one-legged pantyhose with leg extension | ● | ● | ● |
| BT / B'T / CT / DT capri pants | ● | ● | ● |
| FT / ET bermuda pants | ● | ● | ● |

● available in standard sizes 1 – 6 and custom-made ● only available in custom-made

Standard colours


Juzo Expert (CCL 1 – 3 in custom-made)



Juzo Expert (standard sizes and CCL 4), Juzo Expert Silver and Juzo Expert Cotton only in colour Almond.


Trend colours · Batik Collection · Iron-on patches

Juzo Expert (CCL 1 – 3 in custom-made)



See the “Product overview” for the current trend colours and Batik Collection. The Batik Collection is available in all trend and standard colours (except Almond) in tie dye in combination with white or black.

Juzo Expert, Juzo Expert Silver and Juzo Expert Cotton (custom-made)



To select a suitable garment and extra options or for further information, we recommend the brochures “Special-purpose garments” for lymphology and scar therapy.

The strong experts with seam



Juzo® Expert Strong

Reliable and strong products for lymphology


Thanks to its high therapeutic pressure and pleasant massaging effect, Juzo Expert Strong is especially suitable for severe lymphological disorders. The firm and strong fabric provides professional support in the treatment of pronounced lymphoedema and lipoedema.



Closed toes



Open toes



Juzo® Expert Strong Silver

Strong with antibacterial effect


Juzo Expert Strong Silver offers a tailored treatment option for the most severe lymphological disorders. The strong fabric with additional knitted-in silver thread provides the highest possible therapeutic pressure as well as an antibacterial effect that effectively combats bacteria on the skin (e.g. Staphylococcus aureus).



Closed toes



Open toes



Juzo® Expert Strong Cotton

Strong with natural fibres

This fabric combines the reliability and strength of Juzo Expert Strong with the exceptional wearing comfort of cotton. It is perfect for patients with sensitive skin, and people who simply prefer wearing cotton.



Closed toes



Open toes



| Model | Juzo® Expert Strong | Juzo® Expert Strong Silver | Juzo® Expert Strong Cotton |
|--|------------------------------|------------------------------|----------------------------|
| Article number / Compression class | 3051 3052 3053 3054 | 3051 3052 3053 3054 | 3052 3053 |
| AD | ● | ● | ● |
| AD with thin silicone border | ● | ● | ● |
| AD with elastic border | ● | ● | ● |
| AF | ● | ● | ● |
| AF with silicone border | ● | ● | ● |
| AF with decorative silicone border | ● | ● | ● |
| AF with lace border | ● | ● | ● |
| AG | ● | ● | ● |
| AG with silicone border | ● | ● | ● |
| AG with decorative silicone border | ● | ● | ● |
| AG with lace border | ● | ● | ● |
| AG with hip attachment | ● | ● | ● |
| AT | ● | ● | ● |
| AT one-legged pantyhose with leg extension | ● | ● | ● |
| BT/B'T/CT/DT capri pants | ● | ● | ● |
| FT/ET bermuda pants | ● | ● | ● |

● only available in custom-made

Standard colours

Juzo Expert Strong (CCL 1 – 3)




Juzo Expert Strong (CCL 4), Juzo Expert Strong Silver and Juzo Expert Strong Cotton only in colour Almond.

Trend colours · Iron-on patches

Juzo Expert Strong (CCL 1 – 3)



Juzo Expert Strong, Juzo Expert Strong Silver and Juzo Expert Strong Cotton



Compression segments leg



Juzo® Expert and Juzo® Expert Strong

Reliable help – also with silver or cotton

Flat knitted compression segments from Juzo are used especially in scar therapy. They are also helpful in special areas of lymphology.

The Juzo Expert Silver and Juzo Expert Strong Silver with a continuous knitted-in silver thread have a natural antibacterial effect. Silver-containing fabrics can therefore restore the natural balance of inflamed or irritated skin as well as reduce itchiness, thereby promoting better healing.

Juzo Expert Cotton is perfect for patients with sensitive skin, and people who simply prefer wearing cotton. The cotton is also good at transporting away moisture, thereby creating a comfortable skin environment.

Compression foot-toe-portions



Juzo® Expert and Juzo® Expert Strong

Tailored treatment – also with silver or cotton

A foot-toe-portion is especially suitable for oedemas in the forefoot and toe area.

Custom-made Juzo Compression Foot-toe-portions are available in Juzo Expert, Juzo Expert Silver, Juzo Expert Cotton, Juzo Expert Strong and Juzo Expert Strong Silver.



| Model | Juzo® Expert | Juzo® Expert Silver | Juzo® Expert Cotton | Juzo® Expert Strong | Juzo® Expert Strong Silver |
|---|--------------|---------------------|---------------------|---------------------|----------------------------|
| Article number / Compression class | 3021 3022 | 3021 3022 | 3021 3022 | 3051 3052 | 3051 3052 |
| BD Calf segment | ● CCL 2 | ● | ● | ● | ● |
| BF Segment | ● | ● | ● | ● | ● |
| BG Segment | ● | ● | ● | ● | ● |
| DG Thigh segment | ● | ● | ● | ● | ● |
| FG Thigh segment | ● CCL 2 | ● | ● | ● | ● |
| ZE Foot-toe-portion with open or closed toes | ● | ● | ● | ● | ● |

● available in standard sizes 1 – 6 and custom-made ● only available in custom-made

Standard colours

Juzo Expert and Juzo Expert Strong (custom-made)




Juzo Expert (standard sizes), Juzo Expert Silver, Juzo Expert Cotton and Juzo Expert Strong Silver only in colour Almond.

Trend colours · Iron-on patches

Juzo Expert and Juzo Expert Strong (custom-made)



Juzo Expert, Juzo Expert Silver, Juzo Expert Cotton, Juzo Expert Strong and Juzo Expert Strong Silver (custom-made)



Compression stump shrinkers (seamless)














Juzo® Dynamic

Reliable and durable – also with silver

In the period immediately after an amputation, the stump is particularly sensitive and painful. Due to their smoothness Juzo Compression Stump Shrinkers are easy to slip on. The soft and elastic border ensures the stump shrinkers can provide optimal treatment both for stumps with a large circumference and those with adipose, oedemic tissue.

The Juzo Dynamic Silver with a continuous knitted-in silver thread has a natural antibacterial effect.




| Model | Juzo® Dynamic | Juzo® Dynamic Silver |
|--|---------------|----------------------|
| Article number / Compression class | 3511 3512 | 3511 3512 |
|  FG short / medium / long / knee exart. Above-knee stump shrinker with hip attachment (also possible in custom-made) | ● | ● |
|  FG short / medium / long / knee exart. Above-knee stump shrinker without hip attachment | ● | ● |
|  FG short / medium / long / knee exart. Above-knee stump shrinker with silicone border | ● | ● |
|  CF short / medium / long Below-knee stump shrinker | ● | ● |
|  CF short / medium / long Below-knee stump shrinker with silicone border | ● | ● |
|  HD Below-knee stump shrinker up to knee after Pirogoff-amputation | ● | ● |
|  SG Stump shrinker arm | ● | ● |
|  SG Stump shrinker arm with thin silicone border | ● | ● |
|  SH Stump shrinker arm with shoulder and strap | ● | ● |
|  SH Stump shrinker arm with extended shoulder part for being attached to the strap of the bra | ● | ● |
|  SH Stump shrinker arm with extended shoulder part and silicone support for being attached to the strap of the bra | ● | ● |

● available in standard sizes I – VI and custom-made ● only available in custom-made

Standard colours

Juzo Dynamic



Juzo Dynamic Silver



Compression stump shrinkers (with seam)



Juzo® Expert



Custom-made treatment – also with silver or cotton


Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement and is very comfortable to wear.

The Juzo Expert Silver with a continuous knitted-in silver thread has a natural antibacterial effect. Silver-containing fabrics can therefore restore the natural balance of inflamed or irritated skin as well as reduce itchiness.

Juzo Expert Cotton is perfect for patients with sensitive skin, and people who simply prefer wearing cotton.

Samples of patient garments





| Model | Juzo® Expert | Juzo® Expert Silver | Juzo® Expert Cotton |
|---|--------------|---------------------|---------------------|
| Article number / Compression class | 3021 3022 | 3021 3022 | 3021 3022 |
|  FG Above-knee stump shrinker with hip attachment | ● | ● | ● |
|  FG Above-knee stump shrinker without hip attachment | ● | ● | ● |
|  FG Above-knee stump shrinker with silicone border | ● | ● | ● |
|  CF Below-knee stump shrinker | ● | ● | ● |
|  CF Below-knee stump shrinker with silicone border | ● | ● | ● |
|  HD Below-knee stump shrinker up to knee after Pirogoff-amputation | ● | ● | ● |
|  DF/DG Suspension sleeve (without lengthwise seam in the hollow of the knee) with lateral silicone-coated stabilizing bands | 3232 ● DF | | |
|  SG Stump shrinker arm | ● | ● | ● |
|  SG Stump shrinker arm with thin silicone border | ● | ● | ● |
|  SH Stump shrinker arm with shoulder and strap | ● | ● | ● |
|  SH Stump shrinker arm with extended shoulder part for being attached to the strap of the bra | ● | ● | ● |
|  SH Stump shrinker arm with extended shoulder part and silicone support for being attached to the strap of the bra | ● | ● | ● |

● available in standard sizes 1 – 6 and custom-made ● only available in custom-made

Standard colours

Juzo Expert (custom-made)



Juzo Expert Silver, Juzo Expert Cotton and suspension sleeve only in colour Almond.

Trend colours

Juzo Expert (custom-made)



See the "Product overview" for the current trend colours.

To select a suitable garment and extra options or for further information, we recommend the brochures "Special-purpose garments" for lymphology and scar therapy.

Compression sleeves (seamless)



Juzo® Soft

Super soft for feel good moments

The fine microfibre fabric with its optimum moisture transport guarantees a comfortable skin environment and a pleasant, soft feeling on the skin. Juzo Soft combines refined fashion with the ultimate in wearing comfort so you can feel your best.



Juzo® Dynamic

Durability you can count on – also with silver

Juzo Dynamic, with its tight and unobtrusive weave, is visually discreet. It offers a perfect fit, and retains its form even during high levels of physical exertion. This makes it the ideal circular knitted sleeve for demanding requirements.

The Juzo Dynamic Silver with a continuous knitted-in silver thread has a natural antibacterial effect.



Juzo® Dynamic Combined Compression Sleeve

Sleeve with knitted-on hand portion

The Juzo Dynamic is also available with a knitted-on gauntlet with thumb stub in reduced compression. The circular knitted hand portion comes with a special wedge in the thumb area and an elastic border with seamless edge on the hand portion. The small silicone border prevents the sleeve from slipping.



| Model | Juzo® Soft | Juzo® Dynamic | Juzo® Dynamic Silver |
|--|--------------|----------------------|----------------------|
| Article number / Compression class | 2001 2002 | 3511 3512 3513 | 3511 3512 3513 |
| CG / CG max Sleeve | ● | ● CCL 1 – 2 | ● |
| CG / CG max Sleeve with thin silicone border | ● | ● CCL 1 – 2 | ● |
| CH / CH max Sleeve with shoulder and strap | ● | ● CCL 1 – 2 | ● |
| CH Sleeve with extended shoulder part for being attached to the strap of the bra | ● | ● | ● |
| CH Sleeve with extended shoulder part and silicone support for being attached to the strap of the bra | ● | ● | ● |
| AGA Sleeve with knitted-on hand portion (with thumb stub) | | ● CCL 1 – 2 | |
| AGA / AGA max Sleeve with knitted-on hand portion (with thumb stub) and thin silicone border | | ● CCL 1 – 2 | |


● available in standard sizes I – VI, I – VI max and custom-made ● only available in custom-made

Standard colours

Juzo Soft



Juzo Dynamic




Juzo Dynamic Silver



Trend colours · Batik Collection

Juzo Soft (custom-made)



See the “Product overview” for the current trend colours and Batik Collection. The Batik Collection is available in all trend and standard colours (except Almond) in tie dye in combination with white or black.

Compression sleeves (with seam)



Juzo® Expert

Staying mobile – also with silver or cotton

Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement and is very comfortable to wear.

The compression sleeves are also available with incorporated silver (Juzo Expert Silver) or cotton (Juzo Expert Cotton).



Juzo® Expert Strong

Reliable and strong – also with silver

Thanks to its high therapeutic pressure and pleasant massaging effect, Juzo Expert Strong is especially suitable for severe lymphological disorders. The firm and strong fabric provides professional support in the treatment of pronounced lymphoedema and lipoedema.

The Juzo Expert Strong Silver with a continuous knitted-in silver thread has a natural antibacterial effect. Silver-containing fabrics can therefore restore the natural balance of inflamed or irritated skin as well as reduce. That effectively combats bacteria on the skin (e.g. Staphylococcus aureus).

Samples of patient garments





| Model | Juzo® Expert | Juzo® Expert Silver | Juzo® Expert Cotton | Juzo® Expert Strong | Juzo® Expert Strong Silver |
|--|----------------------|----------------------|----------------------|----------------------|----------------------------|
| Article number / Compression class | 3021 3022 3023 | 3021 3022 3023 | 3021 3022 3023 | 3051 3052 3053 | 3051 3052 3053 |
| CG Sleeve | ● | ● | ● | ● | ● |
| CG Sleeve with thin silicone border | ● | ● | ● | ● | ● |
| CG Sleeve with overheight, thin silicone border and seam at the outside of the arm Please indicate: right or left | ● CCL 2 | ● | ● | ● | ● |
| CH Sleeve with shoulder and strap | ● | ● | ● | ● | ● |
| CH Sleeve with extended shoulder part for being attached to the strap of the bra | ● | ● | ● | ● | ● |
| CH Sleeve with extended shoulder part and silicone support for being attached to the strap of the bra | ● | ● | ● | ● | ● |
| AGA / AGO / AGG Sleeve with knitted-on hand portion (with thumb stub) or knitted-on glove (with closed or open fingers) | ● | ● | ● | ● | ● |
| AHAG / AHOG / AHGG Sleeve with knitted-on hand portion (with thumb stub) or knitted-on glove (with closed or open fingers) and shoulder and strap | ● | ● | ● | ● | ● |

Standard colours

● available in standard sizes 2 – 4 and custom made ● only available in custom-made


Juzo Expert and Juzo Expert Strong (custom-made)



Juzo Expert (standard sizes), Juzo Expert Silver, Juzo Expert Cotton and Juzo Expert Strong Silver only in colour Almond.

Trend colours · Batik Collection · Iron-on patches

Juzo Expert (custom-made)




See the “Product overview” for the current trend colours and Batik Collection. The Batik Collection is available in all trend and standard colours (except Almond) in tie dye in combination with white or black.

Juzo Expert Strong



Juzo Expert, Juzo Expert Silver, Juzo Expert Cotton, Juzo Expert Strong and Juzo Expert Strong Silver (custom-made)



To select a suitable garment and extra options or for further information, we recommend the brochures “Special-purpose garments” for lymphology and scar therapy.

Compression gloves (seamless)



Juzo® Classic Seamless



Comfortable thanks to its seamless fabric

Juzo Classic Seamless Gloves are setting new standards in oedema therapy. The innovative seamless fabric is breathable, thereby assuring a comfortable skin environment. Being soft and smooth as well, it guarantees an especially high wearing comfort and easy putting-on and taking-off. The 3D comfort zone in the sensitive area between the thumb and index finger prevents the fabric from creasing and protects the skin from irritation. It is designed not to restrict the thumb movement in any way.

Samples of patient garments





| Model | | Juzo® Classic Seamless | |
|---|---|------------------------|------|
| Article number / Compression class | | 2301 | 2302 |
|  | HDA Gauntlet with thumb stub | ● | |
|  | HFA Glove with open fingers | | ● |

● available in standard sizes 2 – 5

Standard colours



Sesame



Almond



Black pepper



Compression gloves (with seam)



Juzo® Expert

Pleasantly soft for tailored treatment – also with silver or cotton

The anatomy of the hand is very complex, so medical compression garments for the hand need to be very sophisticated. Thanks to its exceptional softness and bi-directional elasticity, Juzo Expert provides the greatest possible freedom of movement and is very comfortable to wear.

The compression gloves are also available with incorporated silver (Juzo Expert Silver) or cotton (Juzo Expert Cotton).



Juzo® Expert Strong

Strong for lymphology – also with silver

Specialized Juzo Expert Strong compression gloves, when used to treat oedemas, prevent lymphatic fluid from accumulating in the hand and promote the return flow of lymph liquid out of the tissues. Thanks to its high therapeutic pressure and pleasant massaging effect, Juzo Expert Strong is especially suitable for severe lymphological disorders.

The Juzo Expert Strong Silver with a continuous knitted-in silver thread has a natural antibacterial effect. Silver-containing fabrics can therefore restore the natural balance of inflamed or irritated skin as well as reduce itchiness.



| Model | Juzo® Expert | Juzo® Expert Silver | Juzo® Expert Cotton | Juzo® Expert Strong | Juzo® Expert Strong Silver |
|--|--------------|---------------------|---------------------|---------------------|----------------------------|
| Article number / Compression class | 3021 3022 | 3021 3022 | 3021 3022 | 3051 3052 | 3051 3052 |
| HDA Short gauntlet (up behind the wrist) with open or closed thumb stub | ● open | ● | ● | ● | ● |
| HDA Long gauntlet (up approx. the middle of the forearm) with open or closed thumb stub | ● | ● | ● | ● | ● |
| HF / HFA Short glove (up behind the wrist) with open or closed fingers | ● | ● | ● | ● | ● |
| HF / HFA Long glove (up approx. the middle of the forearm) with open or closed fingers | ● | ● | ● | ● | ● |

● available in standard sizes 1 – 6 and custom-made ● only available in custom-made

| Model | Juzo® Expert | Juzo® Expert Silver |
|---|--------------|---------------------|
| Article number / Compression class | 3021 3022 | 3021 3022 |
| FI Finger support open or closed | ● open | ● |
| FI Finger support with fastener (velcro) | ● | ● |

● available in standard sizes 1 – 3 and custom-made ● only available in custom-made

Standard colours

Juzo Expert and Juzo Expert Strong (custom-made)



Juzo Expert (standard sizes), Juzo Expert Silver, Juzo Expert Cotton and Juzo Expert Strong Silver only in colour Almond.

Trend colours · Iron-on patches

Juzo Expert and Juzo Expert Strong (custom-made)



See the “Product overview” for the current trend colours.

Juzo Expert, Juzo Expert Silver, Juzo Expert Cotton, Juzo Expert Strong and Juzo Expert Strong Silver (custom-made)



Iron-on patches especially for children’s compression products. Details please see page 9.

To select a suitable garment and extra options or for further information, we recommend the brochures “Special-purpose garments” for lymphology and scar therapy.

Compression segments thorax






Juzo® Expert

Individual treatment – also with silver







Compression thorax segments are used for oedema and scar therapy, e.g. after breast surgery and burns or scalding. These segments can, if necessary, be combined with a pantyhose. Compression bodysuits are particularly popular with children and women.

Custom-made Juzo Compression Thorax Segments are available in Juzo Expert and Juzo Expert Silver.

Samples of patient garments





| Model | Juzo® Expert | Juzo® Expert Silver |
|---|--------------|---------------------|
| Article number / Compression class | 3021 3022 | 3021 3022 |
|  Thorax compression segment from armpit to hip | ● | ● |
|  Thorax compression vest sleeveless slip-on type or with closure option (velcro, zipper or hook / eye fastener) | ● | ● |
|  Thorax compression vest with long sleeves or short sleeves slip-on type or with closure option (velcro, zipper or hook / eye fastener) | ● | ● |
|  Thorax compression bodysuit sleeveless with body fastening only in conjunction with closure option (velcro, zipper or hook / eye fastener) | ● | ● |
|  Thorax compression bodysuit with long sleeves or short sleeves with body fastening only in conjunction with closure option (velcro, zipper or hook / eye fastener) | ● | ● |
|  Gusset sewed-on to one side of the thorax garment | ● | ● |

● only available in custom-made

Standard colours

Juzo Expert




Juzo Expert Silver only in colour Almond.

Trend colours · Iron-on patches

Juzo Expert



Juzo Expert and Juzo Expert Silver



Compression segments head




Juzo® Expert

Aesthetic scar healing – also with silver

Flat knitted fabrics with mild compression 18–21 mmHg (CCL 1), in particular Juzo Expert and Juzo Expert Silver, are recommended for facial compression treatment. Every Juzo compression garment is custom-made based on the measured body dimensions of the wearer. Thus, a fast supply in clinics and hospitals is guaranteed. Juzo offers a variety of facial masks for the neck and head area, and these are tailored according to the location and extend of the scarring.

Custom-made Juzo Compression Segments Head are available in Juzo Expert and Juzo Expert Silver.

Samples of patient garments





| Model | Juzo® Expert | Juzo® Expert Silver |
|---|--------------|---------------------|
| Article number / Compression class | 3021 | 3021 |
| Neck-chin-bandage with head band and neck part closable using velcro or hook fastener at the back of the head | ● | ● |
| Face mask forehead and back of head open openings for ears on request closable using velcro or hook fastener at the back of the head | ● | ● |
| Face mask forehead and back of head closed openings for ears on request closable using velcro or hook fastener at the back of the head | ● | ● |
| Face mask forehead and back of head open openings for eyes, nose, mouth and ears on request closable using velcro or hook fastener at the back of the head | ● | ● |
| Face mask forehead and back of head closed openings for eyes, nose, mouth and ears on request closable using velcro or hook fastener at the back of the head | ● | ● |

● only available in custom-made

Standard colours

Juzo Expert



Juzo Expert Silver only in colour Almond.

Trend colours · Iron-on patches

Juzo Expert



See the "Product overview" for the current trend colours.

Juzo Expert and Juzo Expert Silver



Iron-on patches especially for children's compression products. Details please see page 9.



Anti-thrombosis stockings



Varodem® Hospital

Anti-thrombosis stocking with inspection opening


The Varodem Hospital anti-thrombosis stocking promotes venous return before and after operative treatment. Its main applications are for thrombosis and embolism prophylaxis in immobilized patients confined to bed.

The Varodem Hospital has an inspection opening. This allows constant monitoring of the blood circulation situation. The stocking is made of latex-free yarns. The colour marking around the toes and also near the silicone border indicates the size of the stocking. Thanks to the 3 different lengths the stockings fits to nearly all leg dimensions.

| Model | | Varodem® Hospital |
|---|--|-------------------|
| Article number | | 5601 |
|  | AG short / normal / long with silicone border | ● |
| | AG short / normal / long with hip attachment | ● |

● available in standard sizes S – XXL

Colour



White

Adjustable compression system



Juzo® Compression Wrap

For the treatment of ulcus curis venosum

The Juzo Compression Wrap is a short-stretch garment that can be employed as an alternative to conventional wrap-around bandages. Thanks to its short-stretch properties, it delivers a very high working pressure and a low resting pressure. The velcro fastener makes it easy to apply, and easy to adjust for fluctuations in body circumference. The special indentations prevent „window oedema“ by ensuring that no gaps can form between the straps. The Juzo Compression Wrap is reversible, with black on one side and beige on the other, to suit the individual patient's preference.

| Model | Art . no. | Juzo Compression Wrap |
|--|--------------------|-----------------------|
|  LF normal Foot wrap incl. Stockinette (underliner) | 6000 LF | ● 2 – 4 |
|  LC normal / long Calf wrap incl. Stockinette (underliner) | 6000 LC | ● 2 – 6 |
|  LC normal max / long max Calf wrap incl. Stockinette (underliner) | 6000 LC | ● 2 – 6 |
|  Stockinette (underliner) | 6081 | ● uni |
|  Velcro fasteners set for foot wrap set for calf wrap | 6000 VS 6000 VC | ● uni ● 1 – 2 |

● available in standard / uni sizes


Colour



Samples of patient garments



Adjustable compression system



Juzo® ACS light

Efficient and effective adjustable compression system

The short-stretch compression wraps can be employed as an efficient supplementary treatment during the decongestion phase. The Juzo ACS light products can be adjusted for the individual wearer by means of velcro fasteners. Thanks to the incorporated lining, the wrap is easy to use, can be applied in next to no time and thus increases patient compliance.

Therapeutic compression pressures up to 40 mmHg are possible for arm and leg.

Samples of patient garments



| Model | Art. no. | Standard sizes | Uni size |
|--|----------|----------------|----------|
|  Foot wrap normal / long | J701LI | ● 1 – 5 | |
|  Calf wrap normal / long available for left or right | J702LI | ● 1 – 5 | |
|  Knee wrap normal | J703LI | ● 1 – 5 | |
|  Thigh wrap normal / long available for left or right | J704LI | ● 1 – 5 | |
|  Hand wrap normal | J705LI | ● 1 – 4 | |
|  Arm wrap normal / long available for left or right | J706LI | ● 1 – 4 | |
|  Hip band | J711 | | ● |
|  Midline band supplementary band for the upper extremities | J712 | ● 1 – 2 | |
|  Velcro fasteners | J713 | | ● |
|  Pads set with 2 thicknesses | J714 | ● 1 – 4 | |
|  Anti-slip band 12 x 115 cm | J715 | | ● |
|  Extender | J716 | | ● |

Colour

Black

Juzo® SoftCompress



Juzo® SoftCompress

Special products for oedema therapy

Juzo SoftCompress Products are employed as an effective treatment aid for compression bandaging in oedema treatment. Using the comfortable bandages and accessories can significantly improve the treatment outcome.

The vertically-stitched foam channels create differences in pressure, thereby promoting lymph drainage, while at the same time preventing soft tissue fibrosis. The Juzo SoftCompress products are especially skin-friendly and breathable. The material is washable up to 60 °C.



Juzo® SoftCompress Bandages

Improve therapy results






The comfortable Juzo SoftCompress Bandages can be worn under short-stretch bandages or medical compression stockings, depending on the product that is chosen. They provide additional compression, cushioning and an ideal pressure distribution. The bandages are used for fast and effective bandaging and make self-bandaging easy for the patient. They are available in uni sizes, standard sizes and as custom-made items. Custom-made bandages are ideal for night time treatment.








Juzo® SoftCompress LymphAid and BandageAid

The complete solution for compression bandaging

Juzo SoftCompress LymphAid and BandageAid offer a complete solution for lymphological compression bandaging. Unlike conventional sets, they also include the innovative Juzo SoftCompress materials which can significantly improve the treatment outcome. LymphAid comes with a ready-to-use bandage which is also suitable for self bandaging. The more economical BandageAid contains compression rolls.

| Model | Art. no. | Uni size | Standard sizes | Custom-made |
|--|----------|----------|----------------|-------------|
|  Bandage lower leg | 6830 | | ● 1 – 5 | ● |
|  Bandage lower leg incl. 2 malleo pads | 6831 | ● | | |
|  Bandage lower leg and thigh incl. 2 malleo pads and 1 knee joint pad | 6820 | | | ● |
|  Bandage lower leg and thigh incl. 2 malleo pads, 1 knee joint pad and 2 compression rolls | 6821 | ● | | |
|  Bandage arm incl. 1 hand pad and 1 elbow joint pad | 6822 | | ● 1 – 6 | ● |
|  Bandage arm incl. 1 hand pad and 1 elbow joint pad | 6823 | ● | | |

| Model | Art. no. | Set |
|--|----------|-----|
|  LymphAid leg small 1 bandage lower leg in uni size, 20 gauze bandages, 10 short-stretch bandages and 2 malleo pads | 6800 KL | ● |
|  LymphAid leg large 1 bandage lower leg in uni size, 2 compression rolls, 20 gauze bandages, 14 short-stretch bandages, 2 malleo pads and 1 knee joint pad | 6800 GR | ● |
|  LymphAid arm 1 bandage arm in uni size, 20 gauze bandages, 7 short-stretch bandages, 1 hand pad and 1 elbow joint pad | 6801 | ● |
|  BandageAid leg small 1 compression roll, 20 gauze bandages, 5 short-stretch bandages and 2 malleo pads | 6810 KL | ● |
|  BandageAid leg large 3 compression rolls, 20 gauze bandages, 9 short-stretch bandages, 2 malleo pads and 1 knee joint pad | 6810 GR | ● |
|  BandageAid arm 1 compression roll, 20 gauze bandages, 7 short-stretch bandages, 1 hand pad and 1 elbow joint pad | 6811 | ● |

Juzo® SoftCompress



Juzo® SoftCompress Pads

Targeted increase in pressure

Juzo SoftCompress Pads are employed as an effective treatment aid for oedema. The vertically-stitched foam channels create differences in pressure, thereby promoting lymph drainage, while at the same time preventing soft tissue fibrosis.



Juzo® SoftCompress Compression Roll

Additional pressure














The thin, very soft Juzo SoftCompress Compression Rolls already produce compression pressure if applied suitably. The compression rolls distribute the pressure of the additionally applied short-stretch bandages, thereby avoiding constrictions.



Juzo® SoftCompress Compression Sheet

Individual fitting

Juzo SoftCompress Compression Sheets are very thin and can be cut to size for use in various body regions.

| Model | Art. no. | Standard / Uni sizes | Custom-made |
|--|--------------|----------------------|-------------|
|  Breast pad | 6837 | ● 1–3 | ● |
|  Breast bandage | 6838 | ● 1–6 | |
|  Genital pad women | 6839 | ● 1–3 | ● |
|  Genital pad men | 6840 | ● 1–3 | ● |
|  Knee joint pad 13 x 17 cm | 6845 | ● | ● |
|  Elbow joint pad 10 x 13 cm | 6842 | ● | ● |
|  Malleo pad 6.5 x 9 cm delivered in pairs | 6843 | ● | ● |
|  Hand pad 8.5 x 24 cm | 6844 | ● | ● |
|  Compression roll | 6835 | ● 1–4 | |
|  Compression sheet | 6836 | ● 1–2 | |
|  Short-stretch bandage available per piece or package (10 pieces) 6 / 8 / 10 / 10 cm x 5 m | 6846 | ● 1–4 | |
|  Rapidelast V – gauze bandage delivered per package with 20 pieces 6 cm x 4 m (for fingers) 4 cm x 4 m (for toes) | 6847 6848 | ● ● | |
|  Tubular bandage 6 / 8 / 10 cm x 20 m | 6849 | ● 1–3 | |



Postoperative care of the chest area



comfortable fabric




functional zone

Juzo® Postoperative Compression Bra Julia

Seamless technology for optimal breast care

The Juzo Postoperative Compression Bra Julia optimally stabilizes and supports the operated breast after surgical procedures. The different compression zones with their controlled compressive effect promote the healing process and feeling of well-being of the patient.

The straps and the zipper can be opened fully at the front. This makes follow-up examinations easier and enables physicians and care staff to provide better wound care. To avoid further irritation of the sensitive skin of the operated breast by pressure points, the Juzo Postoperative Compression Bra Julia has been knitted so as to be completely free of seams.



comfortable fabric





functional zone

Juzo® Postoperative Compression Bra Jana

Seamless and metal-free support of the chest area

The Juzo Postoperative Compression Bra Jana, when used in conjunction with a wound dressing, can be worn directly over the operated breast area in the immediate postoperative period. Thanks to its moderate compression, it provides optimal support and care for the breast.

The bra can be opened easily and fully by means of the velcro fastener at the middle of the chest. This reduces any additional strain on the wound and makes it easier for physicians and care staff to care for the wound. The seamless underband ensures no irritating pressure points or rubbing arise under the breast.

| Model | Art. no. |
|--|----------|
|  <p>Juzo® Postoperative Compression Bra Julia with different compression zones, without seams, with incorporated pockets on both sides for bra inserts, with zipper and hook / eye underlap Sizes: S – XXL</p> | 9056 |
|  <p>Juzo® Postoperative Compression Bra Jana with different compression zones, without seams, metal-free, with velcro fastener Sizes: S – XXXL</p> | 9055 |

Colours

Juzo Postoperative Compression Bra Julia and Juzo Post-OP Postoperative Compression Bra Jana




White



Black

Juzo® ScarPad silicone patch



Juzo® ScarPad Silicone patch for scar therapy







Scars can result from a completely normal wound healing process. Juzo ScarPad is a self-adhesive scar treatment patch for the effective and individualised treatment of all types of scars. The patch is versatile to use, the treatment is pain-free with no side effects and also suitable for children. Juzo ScarPad comprised of 100 % medical silicone and a special elastic material.

It improves the appearance of the scar by making it less visible, softer and more elastic. The patch is reusable – one patch lasts for approx. 4 – 8 weeks.



Juzo® ScarPad Fingercaps Silicone compression finger stall

Juzo ScarPad Fingercaps consist of a sewn-on silicone and a textile compression finger stall. It offers an optimal and individual treatment for fingers and partially amputated fingers. The Juzo ScarPad Fingercaps can be cut to any length.

| Model | Art. no. | Packaging |
|--|----------|-----------|
|  Juzo® ScarPad light small thickness 0.4 mm, silicone patch 5 x 8 cm (can be cut to any size) incl. Juzo ScarPad Cleaning soap for the daily cleaning | 9414 | 2 pieces |
|  Juzo® ScarPad light thickness 0.4 mm, silicone patch 10 x 14 cm (can be cut to any size) incl. Juzo ScarPad Cleaning soap for the daily cleaning | 9410 | 2 pieces |
|  Juzo® ScarPad extra thickness 1.4 mm, silicone patch 10 x 14 cm (can be cut to any size) incl. Juzo ScarPad Cleaning soap for the daily cleaning | 9411 | 1 piece |
|  Juzo® ScarPad strong thickness 2.3 mm, silicone patch 14 x 22 cm (can be cut to any size) incl. Juzo ScarPad Cleaning soap for the daily cleaning | 9412 | 1 piece |
|  Juzo® ScarPad Fingercaps silicone finger stall with compression cap. size S – XL (can be cut to any length) incl. Juzo ScarPad Cleaning soap for the daily cleaning | 9415 | 4 pieces |
|  Juzo® ScarPad Cleaning Soap special cleansing agent for use with Juzo ScarPad products 50 ml content per tube | 9413 | 1 tube |


Juzo® Liner Stockings Silver



Juzo® Liner Stockings Silver

With antibacterial effect, versatile

These soft-fabric stockings are especially comfortable to wear and are also suitable as a compression aide for patients with extremely sensitive skin. The stockings offer permanent skin protection during slipping-on and taking-off and, thanks to the knitted-in silver thread, also have an antibacterial action. Juzo Liner Stockings Silver are ideally suited for incorporating wound dressings or securing special pressure pads.

| Model | Juzo® Liner Stockings Silver |
|--|---|
| Article number | 9640 |
|  AD set with 3 liner stockings | size S/M – circumference “cC” < 50 cm size L/XL – circumference “cC” > 50 cm |

Colour



Juzo® Special Socks



Juzo® Allround

Maximum elasticity in the shaft and leg band area

Juzo Allround is especially versatile and comfortable to wear thanks to its high elasticity. This flexible sock is ideal for use by wearers of dressings, supports or orthoses, or for wide feet, high insteps, foot or limb deformations or swelling.


The air comfort sole and a special pad protect the sensitive areas of the foot from, for example, pressure points or chafing-related wounds. Juzo Allround is therefore especially suitable for diabetics and rheumatics.



Juzo® Sensitiv

Seamless toe and elastic-free leg band

Thanks to the use of a cotton core yarn, Juzo Sensitiv is especially light and gentle to wear – ideal for diabetics and rheumatics. The natural and soft stretchability of the material ensures a smooth, crease-free fit. The seamless toe, elastic-free leg band and extremely stretchable body also prevent constriction. Comfort does not get any better than this.

| Model | Juzo® Allround | Juzo® Sensitiv |
|--|----------------|----------------|
| Article number | 9750 | 9751 |
|  AB | ● II – V | ● I – V |

● available in standard sizes

Colours

Juzo Allround



Light grey



Marine



Black

Juzo Sensitiv



Ivory




Marine



Black


Juzo® Light Line – support stockings



Juzo® Light Line 5070 and Juzo® Light Line 5140

Elegant with effective support




Thanks to its silky shine and elegant transparency, the Light Line range is ideal for the fashion conscious woman. This support stocking shapes the legs while also improving venous blood flow, massaging the connective tissue, and keeping the legs fit. There is no better way to preserve your health and prevent venous complaints.



Juzo® Light Line 5800

Comfort through cotton – for men and women


These support stockings are especially soft, offer noticeable relief to the legs, and also effectively prevent venous complaints. They also activate the leg muscles and promote venous function. This makes them the ideal companion during work, leisure, travel and sport – both for men and women.

| Model | Juzo® Light Line (70 denier) | Juzo® Light Line (140 denier) | Juzo® Light Line with cotton |
|--|---------------------------------------|---------------------------------------|---------------------------------------|
| Article number / Support class | 5070 2 (approx. 10-14 mmHg) | 5140 3 (approx. 15-18 mmHg) | 5800 3 (approx. 15-18 mmHg) |
|  AD | ● | ● | ● |
|  AG with lace border | ● | ● | |
|  AT | ● | ● | |

● available in standard sizes I – VI

Colours

Juzo Light Line 5070 and Juzo Light Line 5140



Juzo Light Line 5800



Juzo® Slip-on and Taking-off Aids



Juzo Easy Fit

Art. no. 1890

For compression stockings

with open or closed toes

Package unit 1 piece

Slip-on and taking-off aid

| Models | Ankle circumference |
|-------------------|---------------------|
| Juzo® Easy Fit | < 26 cm |
| Juzo® Easy Fit XL | > 26 cm |



Juzo Slippie Leg

Art. no. 1901

For compression stockings / pantyhose

with open toes

Package unit 20 pieces

Slip-on aid

Made of strong Tyvek, excellent gliding effect



Juzo Slip-on Aid

Art. no. 1900

For compression stockings / pantyhose

with open toes

Package unit 40 pieces

Slip-on aid

Pliable nylon fabric, smooth surface



Juzo Slippie Arm

Art. no. 1902

For compression sleeves and elbow supports

Package unit 1 piece

Slip-on aid

Tearproof material, excellent gliding effect



Juzo Special Gloves

Art. no. 9301

Package unit 1 pair

High-quality textile glove with rubberised palm

Long-life gloves for making donning compression garments easier

| Size | Glove size |
|------|------------|
| S | 5 |
| M | 7 |
| L | 9 |
| XL | 11 |

Juzo® Arion Slip-on and Taking-off Aids



Juzo Arion Sim-Slide

Art. no. 1705

For compression stockings / pantyhose

with open toes

Package unit 1 piece

Smooth, two-layer material

Simply fasten to the ankle using the elastic strap

Increases the independence and therapy compliance of patients

Reduces the physical exertion required by the patient or nursing staff

| Size | Shoe size (UK) |
|------|-------------------|
| M | 38 – 41 (5 – 7.5) |
| L | 42 – 45 (8 – 10) |
| XL | ≥ 46 (10.5) |



Juzo Arion Easy-Slide Leg + Easy-Off

Art. no. 1707

For compression stockings / pantyhose

with open toes

Package unit 1 set

Smooth, two-layer material

With velcro fastener for easy fixation above the ankle

Juzo Arion Easy-Off is also suitable for closed toe

Increases the independence and therapy compliance of patients

Reduces the physical exertion required by the patient or nursing staff

| Size | Shoe size (UK) |
|------|-------------------|
| S | 34 – 37 (2 – 4.5) |
| M | 38 – 41 (5 – 7.5) |
| L | 42 – 45 (8 – 10) |
| XL | ≥ 46 (10.5) |

Juzo® Arion Slip-on Aids



Juzo Arion Easy-Slide Arm

art. no. 1701

For compression sleeves and elbow supports

Package unit 1 piece

Smooth, two-layer material

Increases the independence and therapy compliance of patients

Reduces the physical exertion required by the patient or nursing staff

| Size | Arm length |
|------|------------|
| M | 55.5 cm |
| L | 62.5 cm |



Juzo Arion Easy-Slide Arm

art. no. 1701

For compression sleeves and stockings / pantyhose

with open toes

Package unit 1 piece



Smooth, two-layer material

Increases the independence and therapy compliance of patients

Reduces the physical exertion required by the patient or nursing staff

| Size | Arm / leg length |
|------|------------------|
| XS | 25.5 cm |
| S | 39.5 cm |



Juzo Arion Easy-Slide Leg

art. no. 1702

For compression stockings / pantyhose

with open toes

Package unit 1 piece

Smooth, two-layer material

Increases the independence and therapy compliance of patients

Reduces the physical exertion required by the patient or nursing staff

| Size | Shoe size (UK) |
|------|----------------|
| S | 34 – 37 |
| M | 38 – 41 |
| L | 42 – 45 |
| XL | ≥ 46 |



Juzo Arion Magnide

art. no. 1703

For compression stockings / pantyhose

with closed toes

Package unit 1 piece

Smooth material, with innovative magnetic closure

Increases the independence and therapy compliance of patients

Reduces the physical exertion required by the patient or nursing staff

| Size | Shoe size (UK) |
|------|----------------|
| M | ≤ 40 |
| L | 41 – 45 |
| XL | ≥ 46 |



Juzo Arion SlideX

art. no. 1704

Extension for Juzo Arion Easy-Slide Leg and

Juzo Arion Magnide

Package unit 1 piece

Smooth material

Especially for large circumferences

Reduces the physical exertion required by the patient

| Size | Thigh circumference |
|------|---------------------|
| M | ≤ 40 cm |
| L | 41 – 60 cm |
| XL | 61 – 100 cm |



Juzo Arion Proth-Pro

art. no. 1706

For compression above-knee prosthesis

with suction socket

Package unit 1 piece

Smooth material

With extender for even easier use

Reduces the physical exertion required by the patient

| Size | Circumference | Length |
|------|---------------|---------|
| XS | ≤ 40 cm | 12 cm |
| S | 41 – 45 cm | 13 cm |
| M | 46 – 58 cm | 13 cm |
| L | 59 – 70 cm | 15.5 cm |
| XL | 71 – 80 cm | 16 cm |
| XXL | ≥ 81 cm | 17 cm |

Juzo® Care and Cleaning Products



Juzo Adhesive Lotion

Art. no. 9800

Available as 60 ml roll-on

Special-purpose adhesive for reliable and easy securing of compression fabrics with a top border



Juzo Special Detergent

Art. no. 9000

Available as 30 ml and 250 ml tube

For medical compression garments and elastic fabrics
In a handy tube, which permits exact quantity control
Highly concentrated and economical to use



Juzo Vital Balsam 7

Art. no. 9830

Available as 30 ml and 100 ml tube

Special skin care for compression therapy
Cools, nourishes and revitalizes irritated skin



Juzo Cooling Spray

Art. no. 9870

Available as 150 ml can

Foam spray invigorates and treats stressed skin
For an intense feeling of freshness
Can also be applied directly to compression and support stockings

only available in the EU



Juzo Lymph Lotion

Art. no. 9860

Available as 500 ml tube

Skin care lotion with 4 % urea
For daily use
In the treatment of lymphoedema and lipoedema and also venous afflictions

Juzo® Pressure and Lymph Pads



Pressure pads by Dr. Rotter

Art. no. 9400, 9401*
9409 (silicone, without cover)
Packaging unit 1 pair

Covered with textile or silver-containing (antibacterial) material
For DIY fitting



Pressure pad for the hand

Art. no. short model 9420, 9421*
long model 9420L, 9421L*
Packaging unit 1 piece

Covered with textile or silver-containing (antibacterial) material
To be inserted between back of the hand and glove



Juzo Lymphpad Line

Art. no. 9430
Packaging unit 1 piece

Treatment aid for lymph- and lipoedema
Synthetic latex, polyamide fleece on the rear
Size 20 x 29.5 cm (can be cut to any size)



Juzo Lymphpad Square

Art. no. 9431
Packaging unit 1 piece

Treatment aid for lymph- and lipoedema
Synthetic latex, polyamide fleece on the rear
Size 20 x 29.5 cm (can be cut to any size)



Size charts

The fine fashionable stockings and sporty allround stockings

Juzo® Attractive 2701/2702/2781/2782 • Juzo® Hostess® 2501/2502/2581/2582 • Juzo® Soft 2001/2002 • Juzo® Soft Rib 2021/2022

| Circumferences in cm | Type | Measuring point | I x-small | II small | III medium | IV large | V x-large | VI xx-large |
|----------------------|--------------------------------|-----------------|-------------|-------------|------------|-----------|-----------|-------------|
| | Wide silicone border* | cG | 49 – 58 | 53 – 63 | 57 – 68 | 64 – 75 | 67 – 79 | 74 – 85 |
| | Silicone border** | cG | 41 – 48 | 45 – 52 | 47 – 56 | 54 – 63 | 57 – 66 | 63 – 73 |
| | Leggins | cG | 41 – 58 | 45 – 63 | 47 – 68 | 54 – 75 | 57 – 79 | 63 – 85 |
| | Pantyhose | cF | 41 – 55 | 45 – 60 | 47 – 65 | 54 – 72 | 57 – 76 | 63 – 82 |
| | Lace border | cE | 30 – 37 | 33 – 41 | 35 – 45 | 39 – 48 | 41 – 51 | 44 – 55 |
| | Comfort silicone border | cD | 27 – 33 | 29 – 36 | 31 – 40 | 34 – 43 | 36 – 46 | 39 – 48 |
| | Top border | cC | 29 – 34 | 32 – 38 | 34 – 43 | 38 – 47 | 41 – 50 | 44 – 53 |
| | | cB' | 23 – 27 | 25 – 30 | 27 – 33 | 30 – 36 | 32 – 38 | 35 – 41 |
| | | cB | 18 – 20 | 20 – 22 | 22 – 24.5 | 24.5 – 27 | 27 – 30 | 30 – 33 |
| | | cY | 26.5 – 30.5 | 29.5 – 32.5 | 31.5 – 35 | 34 – 37.5 | 36.5 – 41 | 39.5 – 45 |
| | cA | 17.5 – 22 | 19.5 – 24 | 21.5 – 26.5 | 24 – 29 | 26.5 – 33 | 29.5 – 37 | |

*also wide decorative silicone border **also decorative silicone border

| Lengths in cm | Type | Leg length | Type | Leg length |
|---|--------------------|------------|---|-------------------------------|
| lA: 15 – 18 cm Shoe size regular 37 – 42 (UK: 4 – 8) Shoe size Soft Rib 38 – 44 (UK: 5 – 9.5) | lZ: 23.5 – 26.5 cm | | AD AD short AD extra short | 40 – 46 34 – 39 31 – 34 |
| | | | AG / AT / BT AG / AT / BT short AG / AT extra short | 72 – 83 66 – 72 58 – 66 |

The intensive therapy stockings and sporty allround stockings



Juzo® Dynamic 3511/3512/3513 • Juzo® Dynamic Cotton 3511/3512 • Juzo® Dynamic Cotton Rib 3521/3522

| Circumferences in cm | Type | Measuring point | I x-small | II small | III medium | IV large | V x-large | VI xx-large |
|----------------------|------------------------------|-----------------|-------------|-------------|------------|-----------|-----------|-------------|
| | Wide silicone border* | cG | 50 – 56 | 54 – 62 | 58 – 66 | 64 – 73 | 68 – 77 | 73 – 82 |
| | Silicone border** | cG | 43 – 49 | 47 – 53 | 49 – 57 | 55 – 63 | 59 – 67 | 64 – 72 |
| | Pantyhose | cG | 43 – 56 | 47 – 62 | 49 – 66 | 55 – 73 | 59 – 77 | 64 – 82 |
| | Lace border | cF | 43 – 53 | 47 – 59 | 49 – 63 | 55 – 70 | 59 – 74 | 64 – 79 |
| | Top border | cE | 28 – 34 | 30 – 37 | 33 – 41 | 35 – 43 | 37 – 44 | 40 – 48 |
| | | cD | 25 – 30 | 27 – 33 | 30 – 36 | 33 – 39 | 35 – 41 | 37 – 43 |
| | | cC | 27 – 34 | 29 – 37 | 33 – 41 | 36 – 45 | 39 – 47 | 42 – 50 |
| | | cB' | 24 – 26 | 25 – 29 | 27 – 33 | 31 – 36 | 34 – 39 | 36 – 41 |
| | | cB | 18 – 20 | 20 – 22 | 22 – 24.5 | 24.5 – 27 | 27 – 30 | 30 – 33 |
| | | cY | 26.5 – 30.5 | 29.5 – 32.5 | 31.5 – 35 | 34 – 37.5 | 36.5 – 41 | 39.5 – 45 |
| | cA | 17.5 – 22 | 19.5 – 24 | 21.5 – 26.5 | 24 – 29 | 26.5 – 33 | 29.5 – 37 | |

*also wide decorative silicone border **also decorative silicone border

| Lengths in cm | Type | Leg length | Type | Leg length |
|---|--------------------|------------|---|-------------------------------|
| lA: 15 – 18 cm Shoe size regular 37 – 42 (UK: 4 – 8) Shoe size Dynamic Cotton Rib 39 – 45 (UK: 5.5 – 10), short foot part 36 – 39 (UK: 3.5 – 5.5) | lZ: 23.5 – 26.5 cm | | AD AD short AD extra short | 40 – 46 34 – 39 31 – 34 |
| | | | AG / AT AG / AT short AG / AT extra short | 72 – 83 66 – 72 58 – 66 |



The intensive therapy stockings

Juzo® Helasana® 4512

| Circumferences in cm | Type | Measuring point | II small | III medium | IV large | V x-large |
|---|----------------|--------------------|----------------|--------------------|-----------|-------------|
| | Top border | cG | 44 – 66 | 48 – 71 | 52 – 76 | 56 – 81 |
| | | cF | 42 – 62 | 46 – 67 | 50 – 72 | 54 – 77 |
| | | cE | 30 – 40 | 33 – 44 | 36 – 48 | 39 – 52 |
| | | cD | 27 – 37 | 30 – 41 | 33 – 45 | 36 – 49 |
| | | cC | 28 – 38 | 31 – 42 | 34 – 46 | 37 – 50 |
| | | cB' | 23 – 30 | 26 – 34 | 29 – 38 | 32 – 42 |
| | | cB | 18 – 21.5 | 21.5 – 25 | 25 – 28.5 | 28.5 – 31.5 |
| | | cY | 26 – 33 | 29 – 36 | 31 – 39 | 33 – 42 |
| | | cA | 17 – 23.5 | 20.5 – 27 | 24 – 30 | 27 – 33.5 |
| | | Lengths in cm | Type | Leg length | Type | Leg length |
| lA: 15 – 17 cm Shoe size 37 – 41 (UK: 4 – 7.5) | AD AD short | 40 – 46 34 – 40 | AG AG short | 72 – 83 66 – 72 | | |

Ulcer Pro

Juzo® Ulcer Pro 3073 · Juzo® Ulcer Pro Add 2571



| Circumferences in cm | Measuring point | 1 x-small | 2 small | 3 medium | 4 large | 5 x-large | 6 xx-large |
|--|-----------------|--------------------|-------------|-------------|-----------|-------------|-------------|
| | cD | 27 – 29.5 | 29.5 – 32.5 | 32.5 – 35.5 | 35.5 – 39 | 38.5 – 42 | 41.5 – 45.5 |
| | cC | 29 – 31.5 | 31.5 – 34.5 | 34.5 – 37.5 | 37.5 – 41 | 40.5 – 44.5 | 43 – 47 |
| | cB' | 23.5 – 25.5 | 25.5 – 28 | 27.5 – 30 | 30 – 33 | 32.5 – 35.5 | 35 – 38.5 |
| | cB | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 | 28 – 30.5 |
| | cY | 28 – 30.5 | 29 – 31.5 | 31 – 34 | 33 – 36 | 35 – 38.5 | 37 – 40.5 |
| | cA | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 | 28 – 30.5 |
| Lengths in cm | Type | Leg length | | | | | |
| Juzo Ulcer Expert lA: 15 – 18 cm Shoe size 37 – 42 (UK: 4 – 8) | AD AD short | 40 – 46 34 – 39 | | | | | |
| Juzo Ulcer Liner lZ: 23.5 – 26.5 cm | | | | | | | |

All circumference measurements should be taken over the wound dressing (if present)



The experts with seam

Juzo® Expert 3021/3022 · Juzo® Expert Cotton 3021/3022

| Circumferences in cm | Measuring point | 1 x-small | 2 small | 3 medium | 4 large | 5 x-large | 6 xx-large |
|---|-----------------|--------------------|-------------|-------------|-------------|-------------|-------------|
| | cD max | 29.5 – 32.5 | 32.5 – 35.5 | 35.5 – 39 | 38.5 – 42 | 41.5 – 45.5 | 44.5 – 48.5 |
| | cD | 27 – 29.5 | 29.5 – 32.5 | 32.5 – 35.5 | 35.5 – 39 | 38.5 – 42 | 41.5 – 45.5 |
| | cC max | 31.5 – 34.5 | 34.5 – 37.5 | 37.5 – 41 | 40.5 – 44.5 | 43 – 47 | 45.5 – 49.5 |
| | cC | 29 – 31.5 | 31.5 – 34.5 | 34.5 – 37.5 | 37.5 – 41 | 40.5 – 44.5 | 43 – 47 |
| | cB' max | 26 – 28.5 | 28.5 – 31 | 30.5 – 33.5 | 33 – 36.5 | 35 – 38 | 37.5 – 41 |
| | cB' | 23.5 – 25.5 | 25.5 – 28 | 27.5 – 30 | 30 – 33 | 32.5 – 35.5 | 35 – 38.5 |
| | cB | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 | 28 – 30.5 |
| | cY | 28 – 30.5 | 29 – 31.5 | 31 – 34 | 33 – 36 | 35 – 38.5 | 37 – 40.5 |
| | cA | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 | 28 – 30.5 |
| | Lengths in cm | Type | Leg length | | | | |
| lA: 15 – 17 cm Shoe size 37 – 41 (UK: 4 – 7.5), long foot part 40 – 44 (UK: 7 – 9.5) | AD AD short | 40 – 46 34 – 39 | | | | | |

Juzo Expert Cotton is not available in MAX sizes



Compression leg segments

Juzo® Expert 3022

| Lengths in cm | | Circumferences in cm | | | | | |
|---------------|--|----------------------|-----------|-------------|-----------|---------|----------|
| FG | | 1 | 2 | 3 | 4 | 5 | 6 |
| | | x-small | small | medium | large | x-large | xx-large |
| | | 44 – 47 | 47 – 50.5 | 50.5 – 54.5 | 54.5 – 58 | 58 – 61 | 61 – 65 |
| | | 32 – 35 | 35 – 38 | 38 – 41 | 41 – 44 | 44 – 47 | 47 – 51 |

| Lengths in cm | | Circumferences in cm | | | | | |
|---------------|--|----------------------|-------------|-------------|-----------|-------------|-------------|
| BD | | 1 | 2 | 3 | 4 | 5 | 6 |
| | | x-small | small | medium | large | x-large | xx-large |
| | | 27 – 29.5 | 29.5 – 32.5 | 32.5 – 35.5 | 35.5 – 39 | 38.5 – 42 | 41.5 – 45.5 |
| | | 29 – 31.5 | 31.5 – 34.5 | 34.5 – 37.5 | 37.5 – 41 | 40.5 – 44.5 | 43 – 47 |
| | | 23.5 – 25.5 | 25.5 – 28 | 27.5 – 30 | 30 – 33 | 32.5 – 35.5 | 35 – 38.5 |
| | | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 | 28 – 30.5 |

Compression stump shrinkers

Juzo® Dynamic 3511/3512



| Lengths in cm | | Circumferences in cm | | | | | |
|---------------|--|----------------------|---------|---------|---------|---------|----------|
| FG | | I | II | III | IV | V | VI |
| | | x-small | small | medium | large | x-large | xx-large |
| | | 42 – 45 | 45 – 49 | 49 – 53 | 53 – 57 | 57 – 61 | 61 – 65 |
| | | 30 – 32 | 32 – 35 | 35 – 38 | 38 – 41 | 41 – 44 | 44 – 47 |

| Lengths in cm | | Circumferences in cm | | | | | |
|---------------|--|----------------------|---------|---------|---------|---------|----------|
| CF | | I | II | III | IV | V | VI |
| | | x-small | small | medium | large | x-large | xx-large |
| | | 40 – 42 | 42 – 45 | 45 – 48 | 48 – 51 | 51 – 54 | 54 – 57 |
| | | 30 – 32 | 32 – 35 | 35 – 38 | 38 – 41 | 41 – 44 | 44 – 47 |
| | | 28 – 30 | 30 – 33 | 33 – 36 | 36 – 39 | 39 – 42 | 42 – 45 |

Suspension sleeves

Juzo® Expert 3232



| Lengths in cm | | Circumferences in cm | | | | | |
|---------------|--|----------------------|---------|---------|---------|---------|----------|
| DF | | 1 | 2 | 3 | 4 | 5 | 6 |
| | | x-small | small | medium | large | x-large | xx-large |
| | | 36 – 38 | 38 – 40 | 40 – 42 | 42 – 44 | 44 – 46 | 46 – 48 |
| | | 29 – 31 | 31 – 33 | 33 – 35 | 35 – 37 | 37 – 39 | 39 – 41 |
| | | 26 – 28 | 28 – 30 | 30 – 32 | 32 – 34 | 34 – 36 | 36 – 38 |



Compression sleeves

Juzo® Soft 2001 / 2002 • Juzo® Dynamic 3511 / 3512

| Lengths in cm | Circumferences in cm | | | | | |
|---------------|----------------------|-----------|-----------|-------------|-----------|-------------|
| | I | II | III | IV | V | VI |
| CG | x-small | small | medium | large | x-large | xx-large |
| | 30 – 33 | 32 – 35 | 34 – 37 | 36 – 39.5 | 38 – 41.5 | 40 – 44 |
| | 24 – 26 | 26 – 28.5 | 28 – 30.5 | 30 – 33 | 32 – 35 | 34 – 37 |
| | 28 – 30.5 | 30 – 33 | 32 – 35 | 34 – 37 | 36 – 39.5 | 38 – 41.5 |
| | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | 25 – 27.5 | 27 – 29.5 | 29 – 31.5 | 30.5 – 33.5 | 32 – 35 | 33.5 – 36.5 |
| | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | 25 – 27.5 | 27 – 29.5 | 29 – 31.5 | 30.5 – 33.5 | 32 – 35 | 33.5 – 36.5 |
| | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | 14.5 – 16 | 16 – 17.5 | 17.5 – 19 | 19 – 20.5 | 20.5 – 22 | 22 – 24 |

Juzo® Dynamic 3511 / 3512

| Lengths in cm | Circumferences in cm | | | | | |
|----------------|----------------------|-----------|-----------|-------------|-----------|-------------|
| | I | II | III | IV | V | VI |
| AGA | x-small | small | medium | large | x-large | xx-large |
| | 30 – 33 | 32 – 35 | 34 – 37 | 36 – 39.5 | 38 – 41.5 | 40 – 44 |
| | 24 – 26 | 26 – 28.5 | 28 – 30.5 | 30 – 33 | 32 – 35 | 34 – 37 |
| | 28 – 30.5 | 30 – 33 | 32 – 35 | 34 – 37 | 36 – 39.5 | 38 – 41.5 |
| | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | 25 – 27.5 | 27 – 29.5 | 29 – 31.5 | 30.5 – 33.5 | 32 – 35 | 33.5 – 36.5 |
| | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | 25 – 27.5 | 27 – 29.5 | 29 – 31.5 | 30.5 – 33.5 | 32 – 35 | 33.5 – 36.5 |
| | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | 14.5 – 16 | 16 – 17.5 | 17.5 – 19 | 19 – 20.5 | 20.5 – 22 | 22 – 24 |
| | 17.5 – 19 | 19 – 20.5 | 20.5 – 22 | 22 – 23.5 | 23.5 – 25 | 25 – 26.5 |



Juzo® Expert 3022

| Lengths in cm | Circumferences in cm | | |
|---|----------------------|-------------|-------------|
| | 2 small | 3 medium | 4 large |
| CG with overheight, silicone border and seam at the outside of the arm | 27.5 – 31 | 31 – 35 | 35 – 39 |
| | 24.5 – 27.5 | 27.5 – 30.5 | 30.5 – 34.5 |
| | 23 – 25.5 | 25.5 – 28 | 28 – 31 |
| | 22 – 24 | 24 – 27 | 27 – 30 |
| | 17.5 – 20 | 20 – 22 | 22 – 24 |
| | 16 – 18 | 18 – 20 | 20 – 22 |

Compression gloves

Juzo® Classic Seamless 2301 / 2302



| Lengths in cm | Circumferences in cm | Circumferences in cm | | | |
|---------------|----------------------|----------------------|---------------|-------------|--------------|
| | | II small | III medium | IV large | V x-large |
| HDA | | 17.5 – 19.5 | 19.5 – 21 | 21 – 23 | 23 – 25.5 |
| HFA | | 14.5 – 16.5 | 16.5 – 18 | 18 – 20 | 20 – 22 |



Juzo® Expert 3021 / 3022


| Lengths in cm | Circumferences in cm | Circumferences in cm | | | | | |
|---------------|----------------------|----------------------|------------|-------------|------------|--------------|---------------|
| | | 1 x-small | 2 small | 3 medium | 4 large | 5 x-large | 6 xx-large |
| HDA | | cX | cX | cX | cX | cX | cX |
| | | 6.1 | 6.5 | 6.9 | 7.4 | 7.8 | 8.2 |
| | | 17.5 | 19 | 20.5 | 22 | 23.5 | 25 |
| | | 14.5 | 16 | 17.5 | 19 | 20.5 | 22 |

| Lengths in cm | Circumferences in cm | Circumferences in cm | | | | | |
|---------------|----------------------|----------------------|------------|-------------|------------|--------------|---------------|
| | | 1 x-small | 2 small | 3 medium | 4 large | 5 x-large | 6 xx-large |
| HFA | | cX | cX | cX | cX | cX | cX |
| | | 6.1 | 6.5 | 6.9 | 7.4 | 7.8 | 8.2 |
| | | 6.2 | 6.7 | 7.1 | 7.6 | 8 | 8.5 |
| | | 6 | 6.3 | 6.7 | 7.1 | 7.5 | 7.9 |
| | | 5.7 | 6 | 6.3 | 6.7 | 7.1 | 7.4 |
| | | 5.3 | 5.5 | 5.9 | 6.3 | 6.7 | 7.0 |
| | | 17.5 | 19 | 20.5 | 22 | 23.5 | 25 |
| | | 14.5 | 16 | 17.5 | 19 | 20.5 | 22 |

| Lengths in cm | Circumferences in cm | Circumferences in cm | | | | | | | | | | | |
|---------------|----------------------|----------------------|-----------|------------|-----------|-------------|-----------|------------|-----------|--------------|-----------|---------------|-----------|
| | | 1 x-small | | 2 small | | 3 medium | | 4 large | | 5 x-large | | 6 xx-large | |
| HFK | | cX | cZ | cX | cZ | cX | cZ | cX | cZ | cX | cZ | cX | cZ |
| | | 6.1 | 5.1 | 6.5 | 5.5 | 6.9 | 5.9 | 7.4 | 6.4 | 7.8 | 6.8 | 8.2 | 7.2 |
| | | 6.2 | 4.6 | 6.7 | 5.1 | 7.1 | 5.5 | 7.6 | 6.0 | 8 | 6.4 | 8.5 | 6.9 |
| | | 6 | 4.6 | 6.3 | 4.9 | 6.7 | 5.3 | 7.1 | 5.7 | 7.5 | 6.1 | 7.9 | 6.5 |
| | | 5.7 | 4.6 | 6 | 4.8 | 6.3 | 5.1 | 6.7 | 5.5 | 7.1 | 5.9 | 7.4 | 6.2 |
| | | 5.3 | 4.5 | 5.5 | 5.5 | 5.9 | 4.7 | 6.3 | 5.1 | 6.7 | 5.5 | 7.0 | 5.8 |
| | | 17.5 | | 19 | | 20.5 | | 22 | | 23.5 | | 25 | |
| | | 14.5 | | 16 | | 17.5 | | 19 | | 20.5 | | 22 | |

Compression finger supports

Juzo® Expert 3021 / 3022



| Lengths in cm | Circumferences in cm | Circumferences in cm | | | | | |
|---------------|----------------------|----------------------|-----------|------------|-----------|-------------|-----------|
| | | 1 x-small | | 2 small | | 3 medium | |
| FI | | cX | cZ | cX | cZ | cX | cZ |
| | | 6.0 | 5.0 | 6.7 | 5.7 | 7.5 | 6.5 |



Anti-thrombosis stockings

Varodem® Hospital 5601

| Circumferences in cm | Measuring point | S Red | M Green | L Yellow | XL Black | XXL Brown | Type | Leg length |
|----------------------|-----------------|----------|------------|-------------|-------------|--------------|------|------------|
| | cG | 45 – 55 | 50 – 60 | 55 – 65 | 60 – 70 | 65 – 75 | | 64 – 74 |
| | cC | 26 – 32 | 32 – 38 | 38 – 44 | 44 – 50 | 46 – 52 | | 74 – 84 |
| | cB | 18 – 21 | 21 – 24 | 24 – 27 | 27 – 30 | 29 – 32 | | 84 – 94 |

Adjustable compression systems

Juzo® Compression Wrap


| Foot wrap 6000 | Measurements in cm | 2 | 3 | 4 | | |
|----------------|--------------------------|---------|---------|---------|--|--|
| | | small | medium | large | | |
| | Circumference cB | 19 – 24 | 24 – 29 | 29 – 34 | | |
| | Circumference cA¹ | 15 – 23 | 20 – 28 | 25 – 33 | | |


| Calf wrap 6000 | Measurements in cm | 2 | 3 | 4 | 5 | 6 |
|----------------|-----------------------------|---------|---------|---------|---------|----------|
| | | small | medium | large | x-large | xx-large |
| | Circumference cC max | 35 – 41 | 40 – 46 | 44 – 51 | 48 – 55 | 53 – 60 |
| | Circumference cC | 29 – 35 | 34 – 40 | 37 – 44 | 41 – 48 | 46 – 53 |
| | Circumference cB | 18 – 21 | 21 – 24 | 24 – 27 | 27 – 31 | 31 – 35 |
| | Length AD | | | | | |
| | normal | < 40 | < 40 | < 40 | < 40 | < 40 |
| long | > 40 | > 40 | > 40 | > 40 | > 40 | |


Juzo® ACS light

| Foot wrap J701LI | Measurements in cm | 1 | 2 | 3 | 4 | 5 |
|------------------|---|---------|---------|---------|---------|---------|
| | | x-small | small | medium | large | x-large |
| | Circumference cA (over forefoot) | 22 – 24 | 25 – 27 | 28 – 30 | 31 – 34 | 35 – 40 |
| | Length A | | | | | |
| | normal | 16.5 | 18 | 19.5 | 21 | 23 |
| long | 18.5 | 20 | 21.5 | 23 | 25 | |

| Calf wrap J702LI | Measurements in cm | 1 | 2 | 3 | 4 | 5 |
|------------------|---|---------|---------|---------|---------|---------|
| | | x-small | small | medium | large | x-large |
| | Circumference cC (over mid calf) | 35 – 45 | 40 – 52 | 48 – 60 | 55 – 65 | 60 – 70 |
| | Circumference cB (over ankle) | 20 – 27 | 25 – 32 | 30-38 | 37-44 | 45-52 |
| | Length AD | | | | | |
| | normal | 37 | 37 | 42 | 42 | 44 |
| | long | 42 | 42 | 47 | 47 | 49 |

| | | | | | | |
|--|--|----------|----------|----------|----------|----------|
| Knee wrap J703LI  | Measurements in cm | 1 | 2 | 3 | 4 | 5 |
| | | x-small | small | medium | large | x-large |
| | Circumference cE (around the knee) | 40 – 55 | 50 – 65 | 60 – 75 | 70 – 85 | 85 – 100 |
| | Circumference cD (over top of calf) | 35 – 45 | 40 – 52 | 48 – 60 | 55 – 65 | 60 – 70 |

| | | | | | | |
|---|--|----------|-----------|----------|----------|----------|
| Thigh wrap J704LI  | Measurements in cm | 1 | II | 3 | 4 | 5 |
| | | x-small | small | medium | large | x-large |
| | Circumference cG (largest point) | 50 – 65 | 60 – 75 | 70 – 80 | 75 – 90 | 85 – 100 |
| | Circumference cF (10 cms above E) | 45 – 52 | 50 – 59 | 60 – 66 | 67 – 73 | 74 – 80 |
| | Length FG | | | | | |
| | normal | 23 | 23 | 23 | 23 | 23 |
| | long | 29 | 29 | 29 | 29 | 29 |

| | | | | | |
|--|--|----------|----------|----------|----------|
| Arm wrap J706LI  | Measurements in cm | 1 | 2 | 3 | 4 |
| | | x-small | small | medium | large |
| | Circumference cG (over upper arm) | 25 – 30 | 31 – 35 | 36 – 41 | 42 – 47 |
| | Circumference cE (over elbow) | 23 – 27 | 28 – 32 | 33 – 37 | 38 – 43 |
| | Circumference cC (over wrist) | 14 – 16 | 17 – 19 | 20 – 22 | 23 – 25 |
| | Length CG | | | | |
| | normal | 42 – 47 | 42 – 47 | 42 – 47 | 42 – 47 |
| | long | 46 – 50 | 46 – 50 | 46 – 50 | 46 – 50 |

| | | | | | |
|--|---|----------|----------|----------|----------|
| Hand wrap J705LI  | Measurements in cm | 1 | 2 | 3 | 4 |
| | | x-small | small | medium | large |
| | Circumference cA (over knuckles) | < 18 | 18 – 20 | 20 – 22 | > 22 |

Juzo® ACS light Midline band J712

| Size | Product measurements |
|------|----------------------|
| 1 | 12 x 125 cm |
| 2 | 20 x 125 cm |

Juzo® ACS light Pads J714

| Size | Product measurements |
|------|----------------------|
| 1 | 9 x 9 cm |
| 2 | 10 x 10 cm |
| 3 | 11 x 11 cm |
| 4 | 12 x 12 cm |

Special products for oedema therapy

Juzo® SoftCompress Bandage Lower Leg 6830

| Lengths in cm | Circumferences in cm | 1 x-small | 2 small | 3 medium | 4 large | 5 x-large |
|---------------|----------------------|-------------|-------------|-------------|-----------|-------------|
| | cD | 27 – 29.5 | 29.5 – 32.5 | 32.5 – 35.5 | 35.5 – 39 | 38.5 – 42 |
| | cC | 29 – 31.5 | 31.5 – 34.5 | 34.5 – 37.5 | 37.5 – 41 | 40.5 – 44.5 |
| | cB' | 23.5 – 25.5 | 25.5 – 28 | 27.5 – 30 | 30 – 33 | 32.5 – 35.5 |
| | cB | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 |
| | cY | 28 – 30.5 | 29 – 31.5 | 31 – 34 | 33 – 36 | 35 – 38.5 |
| | cA | 18.5 – 20 | 20 – 22 | 22 – 24 | 24 – 26 | 26 – 28.5 |

Juzo® SoftCompress Bandage Arm 6822

| Lengths in cm | Circumferences in cm | 1 x-small | 2 small | 3 medium | 4 large | 5 x-large | 6 xx-large |
|---------------|----------------------|-----------|-----------|-----------|-------------|-----------|-------------|
| | cG | 30 – 33 | 32 – 35 | 34 – 37 | 36 – 39.5 | 38 – 41.5 | 40 – 44 |
| | cF | 28 – 30.5 | 30 – 33 | 32 – 35 | 34 – 37 | 36 – 39.5 | 38 – 41.5 |
| | cE | 25 – 27.5 | 27 – 29.5 | 29 – 31.5 | 30.5 – 33.5 | 32 – 35 | 33.5 – 36.5 |
| | cD | 23 – 25 | 24.5 – 27 | 26 – 28.5 | 27.5 – 30 | 29 – 31.5 | 30.5 – 33.5 |
| | cC | 14.5 – 16 | 16 – 17.5 | 17.5 – 19 | 19 – 20.5 | 20.5 – 22 | 22 – 24 |

Juzo® SoftCompress Breast Pad 6837

| Size | Cup |
|------|-----------|
| 1 | Cup A – B |
| 2 | Cup C – D |
| 3 | Cup E – F |

Juzo® SoftCompress Breast Bandage 6838

| Size | Cup |
|------|-------|
| 1 | Cup A |
| 2 | Cup B |
| 3 | Cup C |
| 4 | Cup D |
| 5 | Cup E |
| 6 | Cup F |

Juzo® SoftCompress Genital Pad Women 6839

| Size | Product measurements |
|------|----------------------|
| 1 | 20 x 20 cm (wxh) |
| 2 | 30 x 28 cm (wxh) |
| 3 | 40 x 32 cm (wxh) |

Juzo® SoftCompress Genital Pad Men 6840

| Size | Product measurements |
|------|----------------------|
| 1 | 28 x 20 cm (wxh) |
| 2 | 30 x 23 cm (wxh) |
| 3 | 32 x 26 cm (wxh) |

Juzo® SoftCompress Compression Roll 6835

| Size | Product measurements |
|------|---------------------------|
| 1 | 5 mm, 8 cm x 4 m (wxhxl) |
| 2 | 5 mm, 10 cm x 4 m (wxhxl) |
| 3 | 5 mm, 12 cm x 4 m (wxhxl) |
| 3 | 5 mm, 16 cm x 4 m (wxhxl) |

Juzo® SoftCompress Compression Sheet 6836

| Size | Product measurements |
|------|----------------------|
| 1 | 65 x 65 cm (wxh) |
| 2 | 125 x 100 cm (wxh) |

Postoperative compression bras

Juzo® Postoperative Compression Bra Julia 9056

| Size | Underbreast circumference | | | | | | | | | | | |
|------|---------------------------|----|----|----|----|----|----|----|----|-----|-----|-----|
| | cm | 60 | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 |
| | inch | 24 | 26 | 28 | 30 | 32 | 34 | 36 | 38 | 40 | 42 | 44 |
| AA | | | | | | | | | | | | |
| A | | S | | M | | | L | | XL | | XXL | |
| B | | | | | | | | | | | | |
| C | | | | | | | | | | | | |
| D | | | | | | | | | | | | |
| E | | | | | | | | | | | | |

Juzo® Postoperative Compression Bra Julia 9056

| Size | Underbreast circumference | | | | | | | | | | | | | |
|------|---------------------------|----|----|----|----|----|----|----|-----|-----|-----|-----|------|-----|
| | cm | 65 | 70 | 75 | 80 | 85 | 90 | 95 | 100 | 105 | 110 | 115 | 120 | 125 |
| | inch | 26 | 28 | 30 | 32 | 34 | 36 | 38 | 40 | 42 | 44 | 46 | 48 | 50 |
| AA | | | | | | | | | | | | | | |
| A | S | | M | | L | | XL | | XXL | | | | XXXL | |
| B | | | | | | | | | | | | | | |
| C | | | | | | | | | | | | | | |
| D | | | | | | | | | | | | | | |
| E | | | | | | | | | | | | | | |

Silicone finger stalls

Juzo® ScarPad Fingercaps 9415

| Size | Finger circumference |
|------|----------------------|
| S | 5.5 cm |
| M | 6.5 cm |
| L | 7.5 cm |
| XL | 8.5 cm |

Special socks

Juzo® Allround 9750

| Size | Shoe size (UK) |
|------|---------------------|
| I | – |
| II | 35 – 38 (3 – 5) |
| III | 39 – 42 (6 – 8) |
| IV | 43 – 45 (8.5 – 11) |
| V | 46 – 48 (11.5 – 13) |

Juzo® Sensitiv 9751

| Size | Shoe size (UK) |
|------|---------------------|
| I | 36 – 38 (3.5 – 5) |
| II | 39 – 40 (6 – 6.5) |
| III | 41 – 42 (7 – 8) |
| IV | 43 – 45 (8.5 – 11) |
| V | 46 – 48 (11.5 – 13) |

Support stockings

Juzo® Light Line 5070 and Juzo® Light Line 5140

| Size Light Line | Size | Shoe size (UK) | Weight (kg) | Body type |
|-----------------|------------------|---------------------|-------------|-------------|
| I | XS 34 – 36 | 36 – 38 (3.5 – 5) | 43 – 52 | small |
| II | S 36 – 38 | 37 – 39 (6.5 – 6) | 50 – 60 | medium |
| III | M 40 – 42 | 39 – 40 (6 – 6.5) | 57 – 66 | medium |
| IV | L 42 – 44 | 40 – 42 (6.5 – 8) | 63 – 75 | large |
| V | XL 44 – 46 | 41 – 43 (7.5 – 8.5) | 73 – 90 | extra large |
| VI | XXL 48 – 50 / 52 | 41 – 43 (7.5 – 8.5) | > 80 | extra large |

Juzo® Light Line 5800

| Size Light Line | Shoe size (UK) |
|-----------------|----------------------|
| I | 36 – 38 (3.5 – 5) |
| II | 38 – 40 (5 – 6.5) |
| III | 40 – 42 (6.5 – 8) |
| IV | 42 – 44 (8 – 9.5) |
| V | 44 – 46 (9.5 – 11.5) |
| VI | 46 – 48 (11.5 – 13) |



Size charts

Accessories

Special products

Upper extremities

Lower extremities

Compression therapy

Julius Zorn GmbH

Juliusplatz 1
86551 Aichach
Germany
info@juzo.com

www.juzo.com

