

Relax Order Form

TO ORDER:

<https://eshop.jobst-usa.com>

Email: HMS-Elvarex-Orders@essity.com

Phone: (+1) 800-537-1063

Fax: (+1) 800-835-4325



Patient Name / BSN File # _____ DOB _____

Address _____ Gender M ☐ F ☐

City/State/Zip _____

Diagnosis _____

Doctor / Address _____

City/State/Zip _____

PO#	
Original Order ☐	Reorder w Changes ☐
Exact Reorder ☐	Schema # _____

Fitter Name _____ Fitter # _____ Fitter Phone _____

Fitter Facility _____ Fitter email _____

Ship To Acct # _____ Acct Name _____

Address _____ City _____ State _____ Zip _____

Email* _____ Phone _____ FAX _____

*By choosing communication via email (above), I acknowledge that Personal Health Information associated with this purchase may be transmitted from BSN in a non-encrypted manner.

Bill To Acct # _____ Acct Name _____

Address _____ City _____ State _____ Zip _____

Email _____ Phone _____ FAX _____

☐ Last 4 digits of credit card on file _____ OR ☐ New card - call to provide credit card #

Name on CC _____ Exp. _____ Billing Zip Code _____

Armsleeves

Quantity/Class	CCL 1 (15-20 mmHg*)
Left	
Right	

Style

- ☐ CG Armsleeve
☐ AG Armsleeve w/gauntlet

Color

- ☐ Beige
☐ Rose

Lower Extremities



Quantity/Class	CCL 1 (15-20 mmHg*)	CCL 2 (20-30 mmHg*)
Left (AD and AG)		
Right (AD and AG)		

Style

- ☐ Knee High
☐ Thigh High

Color

- ☐ Beige
☐ Rose



* Design Pressure
 ** If cB is <20cm, cB should be measured with 0 tension
 *** If needed, slightly more tension can be used
 † Add 1 cm to IAB and IAC

Arion Easy-Slide

- The user-friendly donning aid for open toe compression stockings and tights
- The lightweight, smooth material provides ease of donning/application

Size	Shoe size	BNR	UOM / Box	Order Qty
X-Small	≤ 2	7965803	1	
Small	2.5 - 5.5	7965804	1	
Medium	6 - 8	7965802	1	
Large	8.5 - 11	7965902	1	
X-Large	≥ 11.5	7966001	1	



Arion Magnide®

- The ultimate donning aid for closed toe compression stockings and tights
- Unique magnetic closure for effective donning
- Has flexible sizing to fit different limb sizes

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	≤ 7.5	7957100	1	
Large	8 - 12	7957101	1	
X-Large*	≥ 13	7957102	1	




Arion Magnide® 2in1

- New and improved application aid
- For putting on and taking off closed toe compression stockings and tights
- All-in-one product
- User-friendly donning mechanism using own body weight

Size	Shoe size	BNR	UOM / Box	Order Qty
Medium	< 7.5	7957110	1	
Large	≥ 8 - 10	7957111	1	
X-Large*	≥ 10.5	7957112	1	


* XL also suitable for extra wide feet with cA up to 28 cm



Arion Sim-Slide®

- Multi-functional application aid
- For putting on and taking off open toe compression stockings and tights
- All-in-one product

Size	Shoe size	BNR	UOM / Box	Order Qty
Small	< 3 - 6	7063301	1	
Medium	6.5 - 8	7063302	1	
Large	8.5 - 10.5	7063303	1	
X-Large	≥ 11	7063304	1	



Arion Easy-Slide Kids

- Specifically designed to alleviate the problems experienced when donning arm and leg sleeves for children

Size	Age	Length	BNR	UOM / Box	Order Qty
X-Small	0-3 yrs	11.6" (29.5cm)	7957108	1	
Small	3-6 yrs	15.5" (39.5cm)	7957109	1	



Arion SlideX®

SlideX is especially helpful for patients with:

- Highly sensitive skin on the lower legs
- Open wounds with a plaster or bandage
- An above-average sized lower leg, caused by lipedema or lymphedema
- To be used in combination with Easy-Slide, Sim-Slide® and Magnide®

Size	Circumference of lower leg	BNR	UOM / Box	Order Qty
Medium	≤ 15.7" (40cm)	7957106	1	
Large	≤ 23.6" (60cm)	7957107	1	
Extra large	≤ 39.5" (100cm)	7957107	1	



Dycem® Non-Slip Material*

- Improves grip
- Provides stability
- Helps with easier donning of a JOBST® garment

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
One size	7763600	1	



Dycem® is a registered trademark of Dycem Ltd.

JOBST® Donning Gloves*

- Help to create friction and provide grip
- Help to evenly distribute the fabric on the leg

* Caution: This product contains natural rubber latex which may cause allergic reactions.

Size	BNR	UOM / Box	Order Qty
X- Small	7740027	1 pair	
Small	7740028	1 pair	
Medium	7740029	1 pair	
Large	7740030	1 pair	
X-Large	7740031	1 pair	



JOBST® Stocking Donner

- Easy-to-use device to help apply JOBST® Hosiery
- Ideal for people who have limited dexterity or have difficulty reaching their feet
- Use with closed and open toe styles
- Available in 2 sizes: Standard and Large calf

Size	BNR	UOM / Box	Order Qty
Standard	110913	1 donner	
Large	110910	1 donner	

